

Hep G2 whole cell lysate ab166833

2 References

製品の概要

製品名	Hep G2 whole cell lysate
特記事項	<p>Cell line: HepG2 (Human hepatocellular carcinoma).</p> <p><i>Hep G2 whole-cell extract was prepared from cell cultures of the human hepatocellular liver carcinoma HepG2 cell line. These cells are epithelial in morphology and can be induced to exhibit polarization under appropriate culture conditions.</i></p> <p>HepG2 cells are most commonly used for studies of hepatocyte function and this HepG2 whole cell lysate is specifically recommended for studies related to:</p> <ul style="list-style-type: none"> - human liver diseases, cancers and drug therapies - HBV etiology - signaling pathways - apoptosis <p>Because of their ability to differentiate into polarized epithelium, these cells are routinely utilized to investigate intracellular protein trafficking, particularly in relation to human liver diseases. HepG2 is the most commonly used cell line in studies pertaining to the regulation of hepatic protein synthesis, particularly the synthesis of acute phase proteins after the onset of a systemic inflammatory response.</p> <p>In addition, these cells are important for studies related to liver metabolism, liver cancers, liver regeneration, liver cytotoxicity, apoptosis and as a model for hepatitis B virus (HBV) viral etiology.</p>
アプリケーション	適用あり: WB

製品の特性

Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
バッファー	<p>pH: 7.90</p> <p>Constituents: 0.009% PMSF, 0.01% Magnesium chloride, 0.008% DTT, 0.48% HEPES, 20% Glycerol, 0.75% Potassium chloride</p>
ライゼート備考	<p>HepG2 whole cell lysate was prepared from cell cultures of the human hepatocellular liver carcinoma HepG2 cell line. These cells are epithelial in morphology and can be induced to exhibit polarization under</p>

appropriate culture conditions.

背景

Hepatitis B virus (HBV) is the major cause of acute and chronic hepatitis, leading to progressive development of necroinflammatory changes in the liver, which can result in cirrhosis and hepatocellular carcinoma. Although the development of an effective vaccine to prevent HBV infection has shown promising results and should lead to its eventual eradication, antiviral chemotherapy remains the only effective method to prevent the progression of the disease in chronic carriers. Therefore, the development of new antiretroviral agents active against HBV is needed. HepG2 cells have an epithelial morphology and are thought to be a very useful model to study HBV virus replication via transfection. Cells are also used for cancer and apoptosis studies (in particular signalling pathway studies).

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab166833の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors