

### IGFBP6 overexpression 293T lysate (whole cell) ab94189

画像数 2

医薬用外劇物

#### 製品の概要

製品名	IGFBP6 overexpression 293T lysate (whole cell)
特記事項	ab94189 is a 293T cell transfected lysate in which Human IGFBP6 has been transiently over-expressed using a pCMV-IGFBP6 plasmid. The lysate is provided in 1X Sample Buffer.
アプリケーション	適用あり: WB


#### 製品の特性

Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors. Similarity: Contains 1 IGFBP N-terminal domain. Contains 1 thyroglobulin type-1 domain. PTM: O-linked glycans consist of hexose (probably Gal), N-acetylhexosamine (probably GalNAc) and sialic acid residues. O-glycosylated with core 1 or possibly core 8 glycans. O-glycosylated on one site only in the region AA 143-168 in cerebrospinal fluid.

#### アプリケーション


**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab94189の使用に適用されず  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.



SDS-PAGE - IGFBP6 overexpression 293T lysate  
(whole cell) (ab94189)

ab94189 at 15 µg/lane on an SDS-PAGE gel



Western blot - IGFBP6 overexpression 293T lysate  
(whole cell) (ab94189)

**All lanes :** Anti-IGFBP6 antibody (**ab60718**) at 1/500 dilution

**Lane 1 :** IGFBP6 overexpression 293T lysate (whole cell)  
(ab94189)

**Lane 2 :** 293T non-transfected Lysate

Lysates/proteins at 25 µg per lane.

#### Secondary

**All lanes :** Goat Anti-mouse IgG (H and L) HRP conjugated at  
1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors