

### MTMR1 overexpression 293T lysate (whole cell) ab94170

画像数 1

医薬用外劇物

#### 製品の概要

製品名	MTMR1 overexpression 293T lysate (whole cell)
特記事項	ab94170 is a 293T cell transfected lysate in which Human MTMR1 has been transiently over-expressed using a pCMV-MTMR1 plasmid. The lysate is provided in 1X Sample Buffer.
アプリケーション	適用あり: WB

#### 製品の特性


Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	MTMR1 is a member of the myotubularin related family of proteins. Members of this family contain the consensus sequence for the active site of protein tyrosine phosphatases. The MTMR1 gene shares structural similarities with the MTM1 gene, and is also located in the same chromosomal region. Myotubularin related proteins are a large subfamily of protein tyrosine phosphatases (PTPs) that dephosphorylate D3 phosphorylated inositol lipids. MTMR1 is a lipid phosphatase that acts on phosphatidylinositol 3-phosphate and phosphatidylinositol (3,5)-bisphosphate.

#### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab94170の使用に適用されます**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.

#### 画像



ab94170 at 15µg/lane on an SDS-PAGE gel.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors