

Unblocked Protein A Beads ab169993

1 References

製品の概要

製品名	Unblocked Protein A Beads
種交差性	交差種: Mouse, Rabbit, Guinea pig, Dog, Pig
製品の概要	<p>Abcam's Unblocked protein A beads are agarose beads that have been extensively validated in chromatin immunoprecipitation assays (ChIP).</p> <p>These beads are intended for isolation of immune complexes (chromatin and specific antibody) in ChIP experiments. The beads are suitable for immunoprecipitation of rabbit polyclonal Abs, mouse IgG2a, IgG2b and IgA, guinea pig IgG, dog IgG, pig IgG.</p> <p>The beads should be washed before use.</p> <p>Beads supplied as a suspension 1:4. Keep the beads in Liquid suspension during storage as drying will result in reduced performance.</p>

アプリケーション	適用あり: IP, ChIP
----------	----------------

製品の特性

保存方法	Store at +4°C. Please refer to protocols.
------	---

内容	2800 µl
Protein A Beads Unblocked	1 x 2.8ml

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab169993の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IP		Use at an assay dependent concentration.

アプリケーション	Abreviews	特記事項
ChIP		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors