

# Anti-Adenovirus Type 5 antibody ab6982

★★★★★ [9 Abreviews](#) [79 References](#) [画像数 4](#)

### 製品の概要

製品名	Anti-Adenovirus Type 5 antibody
製品の詳細	Rabbit polyclonal to Adenovirus Type 5
由来種	Rabbit
アプリケーション	<b>適用あり:</b> Neutralising, ELISA, WB, ICC/IF
種交差性	<b>交差種:</b> Human adenovirus 5
免疫原	Tissue, cells or virus corresponding to Adenovirus Type 5. Purified Adenovirus type 5
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション


**The Abpromise guarantee** [Abpromise保証は、次のテスト済みアプリケーションにおけるab6982の使用に適用されます](#)  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Neutralising		Use at an assay dependent concentration. This antibody will neutralize Ad5 infection in a bioassay (1/25,000 dilution for 50 % inactivation of Ad5at 1000 particles/ml) - PMID: 14605137
ELISA	★★★★★ (1)	Use at an assay dependent concentration. Recognizes Ad5 immobilized on microtiter plates.
WB	★★★★★ (5)	Use at an assay dependent concentration.
ICC/IF	★★★★★ (1)	Use at an assay dependent concentration.

## ターゲット情報


**関連性** Adenoviruses are DNA viruses that cause upper respiratory tract infections, conjunctivitis, and other infections in humans. Adenoviruses were discovered in 1953. About 47 different types have been identified since then, and about half of them are believed to cause human diseases. Infants and children are most commonly affected by adenoviruses. Adenovirus infections can occur throughout the year, but seem to be most common from fall to spring. Adenoviruses are responsible for 3-5% of acute respiratory infections in children and 2% of respiratory illnesses in adults.

## 画像



ab6982 staining Adenovirus type 5 in human HeLa cells by Immunocytochemistry/ Immunofluorescence. The cells were fixed in paraformaldehyde, permeabilised in 0.1% Triton and then blocked using 0.1% BSA for 16 hours at 4°C. Samples were then incubated with primary antibody at 1/800 for 1 hour at 37°C. The secondary antibody used was a donkey anti-rabbit IgG conjugated to Alexa Fluor® 546 (red) used at a 1/1000 dilution.

Immunocytochemistry/ Immunofluorescence - Anti-Adenovirus Type 5 antibody (ab6982)  
Image courtesy of an anonymous Abreview.




Western blot - Anti-Adenovirus Type 5 antibody (ab6982)

Top blot (A): ab6982 at a 1/4000 dilution staining Adenovirus type 5 capsid proteins by Western blot (ECL). Secondary = goat anti-rabbit (HRP conjugate).

Bottom blot (B): As above, except ab6982 omitted.

Top blot (A): ab6982 at a 1/4000 dilution staining Adenovirus type 5 capsid proteins by Western blot (ECL). Secondary = goat anti-rabbit (HRP conjugate). Bottom blot (B): As above, except ab6982 omitted.



Western blot - Anti-Adenovirus Type 5 antibody (ab6982)

Image courtesy of an anonymous Abreview

**All lanes** : Anti-Adenovirus Type 5 antibody (ab6982) at 1/10000 dilution

**All lanes** : human hepatocytes whole cell lysate

Lysates/proteins at 20 µg per lane.

**Secondary**

**All lanes** : HRP conjugated goat anti-rabbit secondary polyclonal antibody at 1/10000 dilution


Developed using the ECL technique.

Performed under reducing conditions.

**Observed band size:** 130 kDa

**Exposure time:** 20 seconds

Blocked with 5% milk



Western blot - Anti-Adenovirus Type 5 antibody (ab6982)

Image courtesy of an anonymous Abreview.

**All lanes** : Anti-Adenovirus Type 5 antibody (ab6982) at 1/10000 dilution

**Lane 1** : Whole cell lysate prepared from U2OS cells.

**Lane 2** : Whole cell lysate prepared from U2OS cells infected with Adenovirus 5.

Lysates/proteins at 35 µg per lane.

#### Secondary

**All lanes** : Goat anti-rabbit conjugated to IRDye 700DX

Blocking step was performed using 5% milk for 1 hour at 23°C.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors