

# Agarose Anti-6X His tag® antibody ab1231

## 1 References

### 製品の概要

---

製品名	Agarose Anti-6X His tag® antibody
製品の詳細	Agarose Goat polyclonal to 6X His tag®
由来種	Goat
標識	Agarose
アプリケーション	<b>適用あり:</b> IP
種交差性	<b>交差種:</b> Species independent
免疫原	Full length native HHHHHH (6-His)(purified): conjugated to KLH. <a href="#">Run BLAST with</a> <a href="#">Run BLAST with</a>
特記事項	Affinity purified antibodies were coupled to agarose beads using a cyanogen bromide method.  HIS-TAG® is a trademark of EMD Biosciences, Inc.  The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.  If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.1% Sodium azide Constituents: 0.0268% PBS, 0.58% Sodium chloride
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション

---

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab1231の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IP		Use at an assay dependent concentration. Use 20-40 µl of gel slurry per 0.1-1 mg of protein lysate or extract.

## ターゲット情報

---

**関連性**      The H-H-H-H-H-H motif is used as a tag on many recombinant proteins to facilitate purification. His-tags can be fused to the amino- or carboxy- termini of proteins in transfected or transformed cells.

**細胞内局在**      Depends upon the localization of the parent protein tagged with hexahistidine.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors