

Agarose Goat IgG, polyclonal - Isotype Control ab104155

2 References

製品の概要

製品名	Agarose Goat IgG, polyclonal - Isotype Control
標識	Agarose
アプリケーション	適用あり: IP
特記事項	<p>Goat IgG was coupled to agarose beads using a cyanogen bromide method. Goat IgG concentration was determined by extinction coefficient: absorbance at 280 nm of 1.4 equals 1.0mg of IgG, before coupling to agarose beads.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	<p>pH: 6.8</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: 0.58% Sodium chloride, 0.0268% PBS</p>
特記事項(精製)	Goat IgG was coupled to agarose beads using a cyanogen bromide method.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG
細胞内局在	Secreted

アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab104155の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IP		Use at an assay dependent dilution. Use 15-25µl of gel slurry per 0.1 to 1 mg of protein lysate or extract. ab104155 is an Immunoprecipitation (IP) control that is used in conjunction with agarose-immobilize antibodies made in goat.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors