

### Agarose Anti-V5 tag antibody ab1229

★★★★★ [1 Abreviews](#) [17 References](#)

#### 製品の概要

製品名	Agarose Anti-V5 tag antibody
製品の詳細	Agarose Goat polyclonal to V5 tag
由来種	Goat
標識	Agarose
アプリケーション	<b>適用あり:</b> IP
種交差性	<b>交差種:</b> Species independent
免疫原	Full length native GKIPNPLLGLDST (V5) (purified): conjugated to KLH.

 [Run BLAST with](#)

 [Run BLAST with](#)

#### 特記事項

Provided as 400µl of slurry (200µl gel + 200µl buffer). 100µg of the affinity purified antibody is attached to the gel at a ratio of 500µg antibody/ml of gel or 250µg IgG/ml of 50% slurry. Affinity purified antibodies were coupled to agarose beads using a cyanogen bromide method.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.8 Preservative: 0.1% Sodium azide Constituents: 0.0268% PBS, 0.58% Sodium chloride
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル

## アプリケーション

**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab1229の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IP	★★★★★ (1)	Use at an assay dependent concentration. Use at a concentration of 15 - 20 µl of gel slurry per 0.1 to 1 mg of protein lysate or extract.

## ターゲット情報

**関連性** The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids (GKPIPNPLLGLDST), although it has also been used with a shorter 9 amino acid sequence (IPNPLLGLD).

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors