

Anti-APRIL/TNFSF13 antibody ab3681

画像数 2

製品の概要

製品名	Anti-APRIL/TNFSF13 antibody
製品の詳細	Rabbit polyclonal to APRIL/TNFSF13
由来種	Rabbit
アプリケーション	適用あり: IHC-P, WB
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to Human APRIL/TNFSF13 aa 196-210. (Peptide available as ab6094)
ポジティブ・コントロール	HL60 and Jurkat whole cell lysate (ab7899).
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.2 Preservative: 0.02% Sodium azide
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab3681の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 10 µg/ml.
WB		Use a concentration of 2 - 4 µg/ml. Predicted molecular weight: 27 kDa.

ターゲット情報

機能	Cytokine that binds to TNFRSF13B/TACI and to TNFRSF17/BCMA. May be implicated in the regulation of tumor cell growth. May be involved in monocyte/macrophage-mediated immunological processes.
組織特異性	Expressed at high levels in transformed cell lines, cancers of colon, thyroid, lymphoid tissues and specifically expressed in monocytes and macrophages.
配列類似性	Belongs to the tumor necrosis factor family.
翻訳後修飾	The precursor is cleaved by furin.
細胞内局在	Secreted.

画像



Lane 1 : Anti-APRIL/TNFSF13 antibody (ab3681) at 1 µg/ml


Lane 2 : Anti-APRIL/TNFSF13 antibody (ab3681) at 2 µg/ml

Lane 3 : Anti-APRIL/TNFSF13 antibody (ab3681) at 4 µg/ml

All lanes : Jurkat cells

Predicted band size: 27 kDa

Observed band size: 27 kDa



Immunohistochemistry of APRIL in human prostate tissue with APRIL/TNFSF13 antibody at 10 µg/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors