

Anti-ASB9 antibody ab97918

1 References [画像数 3](#)

医薬用外劇物

製品の概要

製品名	Anti-ASB9 antibody
製品の詳細	Rabbit polyclonal to ASB9
由来種	Rabbit
アプリケーション	適用あり: IHC-P, WB, ICC/IF
種交差性	交差種: Human
免疫原	Recombinant protein fragment, corresponding to a region within amino acids 1-253 of Human ASB9 (NP_001026909).
ポジティブ・コントロール	HeLaS3, HepG2 and MOLT4 cell lysates for WB; HeLa cells for IF/ICC
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 1.21% Tris, 0.75% Glycine, 10% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab97918の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
IHC-P		1/100 - 1/1000.
WB		1/500 - 1/3000. Predicted molecular weight: 32 kDa.
ICC/IF		1/100 - 1/200.

ターゲット情報

関連性


ASB9 may be a substrate-recognition component of an SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins.

画像



Immunohistochemical analysis of paraffin-embedded A549 xenograft, using ab97918 antibody at 1:500 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ASB9 antibody (ab97918)



Western blot - Anti-ASB9 antibody (ab97918)

All lanes : Anti-ASB9 antibody (ab97918) at 1/1000 dilution

Lane 1 : HeLaS3 whole cell lysate


Lane 2 : HepG2 whole cell lysate

Lane 3 : MOLT4 whole cell lysate

Lysates/proteins at 30 µg per lane.

Predicted band size: 32 kDa

12% SDS PAGE



Immunocytochemistry/ Immunofluorescence - Anti-ASB9 antibody (ab97918)

Immunofluorescence analysis of paraformaldehyde-fixed HeLa cells, using ab97918 at 1/200 dilution. Lower image merged with DNA probe.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors