

Anti-Bacillus anthracis protective antigen antibody ab21268

1 References

製品の概要

製品名	Anti-Bacillus anthracis protective antigen antibody
製品の詳細	Rabbit polyclonal to Bacillus anthracis protective antigen
由来種	Rabbit
アプリケーション	適用あり: ELISA
種交差性	交差種: Bacillus anthracis
免疫原	Synthetic peptide corresponding to 14 amino acids at the carboxy-terminus of the Anthrax protective antigen protein (Genbank accession no. AAF86457).
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab21268の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		

追加情報

ELISA: Use at a concentration of 1 µg/ml.
 Note: at this concentration ab21268 will detect 10 ng of free peptide.

Not tested in other applications.
 Optimal dilutions/concentrations should be determined by the end user.

ターゲット情報

関連性

Anthrax toxins are composed of three distinct proteins, a protective antigen (PA), a lethal factor (LF) and an edema factor (EF). None of these proteins are toxic by themselves and several studies indicate that the anthrax toxin has the familiar A-B enzymatic binding structure, with PA acting as the binding domain and EF and/or LF acting as the active fragments.

細胞内局在

Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors