

Anti-Biotin antibody, prediluted ab19218

1 References

製品の概要

製品名	Anti-Biotin antibody, prediluted
製品の詳細	Goat polyclonal to Biotin, prediluted
由来種	Goat
特異性	By immunoelectrophoresis and ELISA this antibody reacts specifically with KLH-biotin and BSA-biotin. There is no reaction with free KLH or BSA.
アプリケーション	適用あり: ELISA, WB, ICC
種交差性	交差種: Species independent
免疫原	Chemical/ Small Molecule corresponding to Biotin.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 6.80 Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab19218の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		1/1000 - 1/30000. For coating plates: 1/100 - 1/500.
WB		1/1000 - 1/30000. Predicted molecular weight: 1 kDa.

アプリケーション	Abreviews	特記事項
ICC		1/200 - 1/2000.

ターゲット情報

関連性 Biotin is a water soluble vitamin, generally classified as a B complex vitamin, also called vitamin B4. After the initial discovery of biotin, nearly forty years of research were required to establish it as a vitamin. Biotin is required by all organisms but can only be synthesized by bacteria, yeasts, molds, algae, and some plant species. Biotin is required as prosthetic group of enzymes involved in incorporation of carbon dioxide into organic compounds. Biotin has a MW of 244 Da.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors