

Biotin Anti-Filamentous Phages M13 fd F1 antibody [B62-FE2] ab20337

★★★★★ 1 Abreviews

製品の概要

| | |
|----------|---|
| 製品名 | Biotin Anti-Filamentous Phages M13 fd F1 antibody [B62-FE2] |
| 製品の詳細 | Biotin Mouse monoclonal [B62-FE2] to Filamentous Phages M13 fd F1 |
| 由来種 | Mouse |
| 標識 | Biotin |
| アプリケーション | 適用あり: Dot blot, WB, ELISA |
| 種交差性 | 交差種: Enterobacteria phage M13 |
| 免疫原 | Tissue, cells or virus corresponding to Escherichia coli Filamentous Phages M13 fd F1. fd phages from E. coli F+ strain (JM109) |
| エピトープ | Binds to an epitope on pVIII (phage coat protein) covering the N-terminal region of g8p AEGDDPAKAAFDLSLQASAT. |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|-------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー | Preservative: 0.09% Sodium azide Constituents: PBS, 0.5% BSA |
| 精製度 | Protein A purified |
| ポリ/モノ | モノクローナル |
| クローン名 | B62-FE2 |
| ミエローマ | unknown |

| | |
|--------|---------|
| アイソタイプ | IgG2b |
| 軽鎖の種類 | unknown |

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab20337の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| Dot blot | | Use at an assay dependent concentration. |
| WB | | Use at an assay dependent concentration. Can also be used in Immunoassays for the identification of recombinant antigen or antibody phages. |
| ELISA | ★★★★★ (1) | 1/500 - 1/1000. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors