

# Biotin Rat IgG2b, kappa monoclonal [RTK4530] - Isotype control ab18542

### 製品の概要

製品名	Biotin Rat IgG2b, kappa monoclonal [RTK4530] - Isotype control
標識	Biotin
特記事項 (標識)	The solution is free of unconjugated biotin.
アプリケーション	<b>適用あり:</b> Flow Cyt (Intra), Flow Cyt
免疫原	Chemical/ Small Molecule corresponding to Rat conjugated to keyhole limpet haemocyanin. Trinitrophenol + KLH
特記事項	<p>Isotype control has been screened on a variety of resting, activated, live, and fixed mouse and human tissues.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Do Not Freeze.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	モノクローナル
クローン名	RTK4530
アイソタイプ	IgG2b
軽鎖の種類	kappa

## アプリケーション

**The Abpromise guarantee**      Abpromise保証は、 次のテスト済みアプリケーションにおけるab18542の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt (Intra)		
Flow Cyt		

### 追加情報

Flow Cyt: Use at an assay dependent dilution. This product is also suitable for use in intracellular flow cytometry.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors