

Biotin Anti-RELM alpha/Retnla antibody ab39628

3 References

製品の概要

製品名	Biotin Anti-RELM alpha/Retnla antibody
製品の詳細	Biotin Rabbit polyclonal to RELM alpha/Retnla
由来種	Rabbit
標識	Biotin
アプリケーション	適用あり: Flow Cyt, Sandwich ELISA, WB
種交差性	交差種: Mouse
免疫原	Recombinant full length protein corresponding to Mouse RELM alpha/Retnla.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Lyophilized:Reconstitute to a concentration of 50µg/ml with sterile PBS solution containing 0.1% BSA.
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Constituents: 0.1% BSA, PBS
精製度	Immunogen affinity purified
特記事項 (精製)	Purified by affinity chromatography and then biotinylated.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab39628の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use at an assay dependent concentration. PubMed: 21825018 ab200208 - Rabbit polyclonal IgG, is suitable for use as an isotype control with this antibody.
Sandwich ELISA		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

ターゲット情報

機能	Probable hormone.
組織特異性	Highest levels in adipose tissue.
配列類似性	Belongs to the resistin/FIZZ family.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors