

Biotin Anti-Streptococcus pneumoniae antibody ab69534

1 References

製品の概要

製品名	Biotin Anti-Streptococcus pneumoniae antibody
製品の詳細	Biotin Rabbit polyclonal to Streptococcus pneumoniae
由来種	Rabbit
標識	Biotin
特異性	ab69534 reacts with Streptococcus pneumoniae in bacterial and infected tissue samples.
アプリケーション	適用あり: ELISA, ICC/IF
種交差性	交差種: Streptococcus pneumoniae
免疫原	Tissue, cells or virus corresponding to Streptococcus pneumoniae. Streptococcus whole-cell blend (serotypes: 3, 4, 6, 7, 9, 14, 18, 19 & 23).
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS
精製度	Affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab69534の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.

ターゲット情報

関連性 Streptococcus pneumoniae are Gram positive diplococci , alpha hemolytic aerotolerant anaerobes. S. pneumoniae is a significant human pathogen, being recognized as a major cause of pneumonia in the late 19th century, and is currently the subject of many humoral immunity studies. S. pneumoniae is a transient member of the normal flora, colonizing the nasopharynx of 40% of healthy adults and children with no adverse effects. There are more than 90 serotypes: serotypes 6, 14, 18, 19, and 23 are the most prevalent accounting for 60-80% of infections depending on the area of the world. Pneumococcal infection accounts for more deaths than any other vaccine preventable bacterial disease. S. pneumoniae can express several different virulence factors which contribute to some of the clinical manifestations seen during infection.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors