

# Anti-BOTOX antibody ab20641

★★★★★ [1 Abreviews](#) [5 References](#)

### 製品の概要

製品名	Anti-BOTOX antibody
製品の詳細	Rabbit polyclonal to BOTOX
由来種	Rabbit
特異性	Clostridium botulinum Type A Neurotoxin.
アプリケーション	<b>適用あり:</b> ELISA, WB
種交差性	<b>交差種:</b> Species independent
免疫原	Full length protein corresponding to BOTOX.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 2.68% PBS, 20% Glycerol (glycerin, glycerine)
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab20641の使用に適用されます  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		1/10000.
WB	★★★★★ (1)	Use at an assay dependent dilution. Predicted molecular weight: 149 kDa.

## ターゲット情報

**関連性** Clostridium botulinum neurotoxins are zinc-dependent proteases that block the release of acetylcholine from peripheral cholinergic nerve endings by targeting specific membrane proteins including synaptosomal-associated proteins (SNAP), vesicle-associated membrane proteins (VAMP), and syntaxins. Neurotoxin A is characteristically associated with human botulism and known to target and hydrolyse the 197-Gln-Arg-198 bond of SNAP-25.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors