

### Anti-Calretinin antibody [6B3] ab82198

★★★★★ 1 Abreviews

#### 製品の概要

---

製品名	Anti-Calretinin antibody [6B3]
製品の詳細	Mouse monoclonal [6B3] to Calretinin
由来種	Mouse
特異性	No cross-reaction with Calbindin or Calmodulin.
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Rat, Human
免疫原	Recombinant Human Calretinin 22kDa
ポジティブ・コントロール	Cerebral cortex, pleural mesothelioma.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.30 Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Tissue culture supernatant
ポリ/モノ	モノクローナル
クローン名	6B3
アイソタイプ	IgG

**The Abpromise guarantee** **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab82198の使用に適用されます  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		

**追加情報**

WB: 1/40 - 1/50. Predicted molecular weight: 31 kDa.

The antibody is suitable for the detection of Calretinin in neurons and may be used as an aid in differential diagnosis of pleura-mesothelioma. This marker is frequently expressed in epitheloid mesothelioma, while occasionally it turns positive in sarcomatous mesothelioma and carcinoma. It was found to be negative in sarcoma and in 90% of adenocarcinoma.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**ターゲット情報**

<b>機能</b>	Calretinin is a calcium-binding protein which is abundant in auditory neurons.
<b>組織特異性</b>	Brain.
<b>配列類似性</b>	Belongs to the calbindin family. Contains 6 EF-hand domains.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors