

Anti-CD163 antibody [BLR087G] ab265592

リコンビナント

画像数 4

製品の概要

製品名	Anti-CD163 antibody [BLR087G]
製品の詳細	Rabbit monoclonal [BLR087G] to CD163
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	Synthetic peptide within Human CD163 aa 850-900. The exact sequence is proprietary. NP_004235.4 Database link: Q86VB7
ポジティブ・コントロール	IHC-P: Human liver carcinoma and lung carcinoma tissue. Human lung carcinoma pleural fluid.
特記事項	This product is sold under License from Bethyl Laboratories, Inc.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 7.8 Preservative: 0.09% Sodium azide Constituents: 99% Borate buffered saline, 0.1% BSA
精製度	Immunogen affinity purified
特記事項 (精製)	Purified from cell culture supernatant.
ポリ/モノ	モノクローナル
クローン名	BLR087G
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab265592の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
IHC-P		1/100 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報


機能	<p>Acute phase-regulated receptor involved in clearance and endocytosis of hemoglobin/haptoglobin complexes by macrophages and may thereby protect tissues from free hemoglobin-mediated oxidative damage. May play a role in the uptake and recycling of iron, via endocytosis of hemoglobin/haptoglobin and subsequent breakdown of heme. Binds hemoglobin/haptoglobin complexes in a calcium-dependent and pH-dependent manner. Exhibits a higher affinity for complexes of hemoglobin and multimeric haptoglobin of HP*1F phenotype than for complexes of hemoglobin and dimeric haptoglobin of HP*1S phenotype. Induces a cascade of intracellular signals that involves tyrosine kinase-dependent calcium mobilization, inositol triphosphate production and secretion of IL6 and CSF1. Isoform 3 exhibits the higher capacity for ligand endocytosis and the more pronounced surface expression when expressed in cells. After shedding, the soluble form (sCD163) may play an anti-inflammatory role, and may be a valuable diagnostic parameter for monitoring macrophage activation in inflammatory conditions.</p>
組織特異性	<p>Expressed in monocytes and mature macrophages such as Kupffer cells in the liver, red pulp macrophages in the spleen, cortical macrophages in the thymus, resident bone marrow macrophages and meningeal macrophages of the central nervous system. Expressed also in blood. Isoform 1 is the lowest abundant in the blood. Isoform 2 is the lowest abundant in the liver and the spleen. Isoform 3 is the predominant isoform detected in the blood.</p>
配列類似性	<p>Contains 9 SRCR domains.</p>
ドメイン	<p>The SRCR domain 3 mediates calcium-sensitive interaction with hemoglobin/haptoglobin complexes.</p>
翻訳後修飾	<p>A soluble form (sCD163) is produced by proteolytic shedding which can be induced by lipopolysaccharide, phorbol ester and Fc region of immunoglobulin gamma. This cleavage is dependent on protein kinase C and tyrosine kinases and can be blocked by protease inhibitors. The shedding is inhibited by the tissue inhibitor of metalloproteinase TIMP3, and thus probably induced by membrane-bound metalloproteinases ADAMs. Phosphorylated.</p>
細胞内局在	<p>Secreted and Cell membrane. Isoform 1 and isoform 2 show a lower surface expression when expressed in cells.</p>

画像




Immunohistochemical analysis of formalin-fixed, paraffin-embedded human lung carcinoma pleural fluid tissue labeling CD163 with ab265592 at 1/100 dilution. Secondary: HRP-conjugated goat anti-rabbit IgG. Substrate: DAB.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD163 antibody
[BLR087G] (ab265592)



Immunohistochemical analysis of formalin-fixed, paraffin-embedded human lung carcinoma tissue labeling CD163 with ab265592 at 1/100 dilution. Secondary: HRP-conjugated goat anti-rabbit IgG. Substrate: DAB.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD163 antibody
[BLR087G] (ab265592)



Immunohistochemical analysis of formalin-fixed, paraffin-embedded human liver carcinoma tissue labeling CD163 with ab265592 at 1/100 dilution. Secondary: HRP-conjugated goat anti-rabbit IgG. Substrate: DAB.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD163 antibody
[BLR087G] (ab265592)

Why choose a recombinant antibody?



Anti-CD163 antibody [BLR087G] (ab265592)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors