

Anti-CD3 epsilon antibody [C3e/3125R] ab270716

リコンビナント

画像数 3

製品の概要

製品名	Anti-CD3 epsilon antibody [C3e/3125R]
製品の詳細	Rabbit monoclonal [C3e/3125R] to CD3 epsilon
由来種	Rabbit
アプリケーション	適用あり: IHC-P, WB
種交差性	交差種: Human
免疫原	Recombinant fragment within Human CD3 epsilon aa 23-119. The exact sequence is proprietary. Database link: P07766
ポジティブ・コントロール	IHC-P: Human tonsil tissue. WB: MOLT-4 cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Preservative: 0.05% Sodium azide Constituents: PBS, 0.05% BSA
精製度	Protein A/G purified
特記事項(精製)	Ab Purified from Bioreactor Concentrate.
ポリ/モノ	モノクローナル
クローン名	C3e/3125R
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab270716の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 1 - 2 µg/ml. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 24 kDa.

ターゲット情報


機能	The CD3 complex mediates signal transduction.
配列類似性	Contains 1 Ig-like (immunoglobulin-like) domain. Contains 1 ITAM domain.
細胞内局在	Membrane.

画像



Anti-CD3 epsilon antibody [C3e/3125R] (ab270716) at 2 µg/ml + MOLT-4 cell lysate.


Predicted band size: 24 kDa



Formalin-fixed, paraffin-embedded human tonsil tissue stained for CD3 epsilon using ab270716 at 2 µg/ml in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-CD3 epsilon antibody [C3e/3125R] (ab270716)

Why choose a recombinant antibody?

 <p>Research with confidence Consistent and reproducible results</p>	 <p>Long-term and scalable supply Recombinant technology</p>
 <p>Success from the first experiment Confirmed specificity</p>	 <p>Ethical standards compliant Animal-free production</p>

Anti-CD3 epsilon antibody [C3e/3125R] (ab270716)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors