

### Anti-CD48 antibody [341] - BSA and Azide free ab277141

リコンビナント

画像数 3

#### 製品の概要

製品名	Anti-CD48 antibody [341] - BSA and Azide free
製品の詳細	Rabbit monoclonal [341] to CD48 - BSA and Azide free
由来種	Rabbit
アプリケーション	<b>適用あり:</b> Flow Cyt, ICC
種交差性	<b>交差種:</b> Mouse
免疫原	Recombinant fragment corresponding to Mouse CD48 aa 1-216. Genbank: NP_031675.1 Database link: <a href="#">P18181</a>
ポジティブ・コントロール	Flow cyt: Mouse BABL/c splenocytes. ICC: Mouse Spleen cells.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Constituent: 100% PBS  0.2 µm filtered solution
キャリア・フリー	はい
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	341
アイソタイプ	IgG

#### アプリケーション


**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab277141の使用に適用されます**  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		1/25 - 1/100.
ICC		1/20 - 1/100.


## ターゲット情報

<b>機能</b>	Ligand for CD2. Might facilitate interaction between activated lymphocytes. Probably involved in regulating T-cell activation.
<b>配列類似性</b>	Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Contains 1 Ig-like V-type (immunoglobulin-like) domain.
<b>細胞内局在</b>	Cell membrane.

## 画像



Flow cytometric analysis of CD48 expression on mouse BABL/c splenocytes. Cells were stained with ab277141 at 1/25 dilution, then a FITC-conjugated second step antibody (Black). Cells were also stained with an isotype control (Grey). The fluorescence histograms were derived from gated events with the forward and side light-scatter characteristics of intact cells.



Immunocytochemistry - Anti-CD48 antibody  
(ab277141)

Immunofluorescence staining of CD48 in mouse spleen cells. Cells were fixed with 4% PFA, blocked with 10% serum, and incubated with ab277141 at 1/60 dilution at 37°C for 1 hour. Then cells were stained with the Alexa Fluor® 488-conjugated Goat Anti-rabbit IgG secondary antibody (green).

### Why choose a recombinant antibody?



Anti-CD48 antibody [341] - BSA and Azide free  
(ab277141)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors