



# Anti-Clostridium difficile Toxin A antibody ab93318

### 製品の概要

製品名	Anti-Clostridium difficile Toxin A antibody
製品の詳細	Rabbit polyclonal to Clostridium difficile Toxin A
由来種	Rabbit
アプリケーション	<b>適用あり:</b> ELISA
種交差性	<b>交差種:</b> Clostridium difficile
免疫原	Synthetic peptide: C-TIDGKKYYFN conjugated to KLH, corresponding to amino acids 2351-2360 of Clostridium difficile Toxin A  <a href="#">Run BLAST with</a>  <a href="#">Run BLAST with</a>

### 特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Protein G purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab93318の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.

## ターゲット情報

**関連性** Clostridium difficile, or C. difficile, a gram positive anaerobic bacterium, is the major causative agent of colitis and diarrhea that can occur following antibiotic intake. C. difficile infection is one of the most common hospital acquired infections. When C. difficile colonize the gut, they release two potent toxins, toxin A and toxin B, which bind to certain receptors in the lining of the colon and ultimately cause diarrhoea and inflammation of the large intestine.

**細胞内局在** Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors