


Anti-Collagen XV alpha 1 antibody ab58717

1 References

製品の概要

| | |
|----------|---|
| 製品名 | Anti-Collagen XV alpha 1 antibody |
| 製品の詳細 | Rabbit polyclonal to Collagen XV alpha 1 |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: WB, ELISA |
| 種交差性 | 交差種: Human 交差が予測される動物種: Mouse  |
| 免疫原 | Synthetic peptide corresponding to C terminal residues of human Collagen XV alpha 1. |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|----------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| バッファー | pH: 7.2 Preservative: 0.01% Sodium azide Constituents: PBS, 50% Glycerol (glycerin, glycerine) |
| 精製度 | Immunogen affinity purified |
| 特記事項(精製) | Purified by antigen specific affinity chromatography. |
| ポリモノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab58717の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| WB | | Use at an assay dependent concentration. |
| ELISA | | Use at an assay dependent concentration. |

ターゲット情報

| | |
|-------|---|
| 機能 | Structural protein that stabilizes microvessels and muscle cells, both in heart and in skeletal muscle. Restin potently inhibits angiogenesis. |
| 組織特異性 | Expressed predominantly in internal organs such as adrenal gland, pancreas and kidney. |
| 配列類似性 | Belongs to the multiplexin collagen family. Contains 4 collagen-like domains. Contains 1 laminin G-like domain. |
| 翻訳後修飾 | Prolines at the third position of the tripeptide repeating unit (G-X-Y) are hydroxylated in some or all of the chains. O-glycosylated; with core 1 or possibly core 8 glycans. An N-terminal peptide contains chondroitin sulfate. |
| 細胞内局在 | Secreted, extracellular space, extracellular matrix. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors