

Anti-DARPP32 antibody [EP721Y] ab40802

リコンビナント **RabMAb**

3 References **画像数 6**

製品の概要

製品名	Anti-DARPP32 antibody [EP721Y]
製品の詳細	Rabbit monoclonal [EP721Y] to DARPP32
由来種	Rabbit
アプリケーション	適用あり: ICC/IF, WB, IHC-P, IP
種交差性	交差種: Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	Human brain lysate and human gastric adenocarcinoma. SH-SY5Y (human neuroblastoma) and human fetal forebrain tissue lysates.
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p> <p>Mouse, Rat: We have preliminary internal testing data to indicate this antibody may not react with these species. Please contact us for more information.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.20</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 40% Glycerol (glycerin, glycerine), 0.11% BSA, 59% PBS</p>
精製度	Protein A purified
ポリ/モノ	モノクローナル

クローン名	EP721Y
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab40802の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		1/50 - 1/100.
WB		1/5000 - 1/50000. Detects a band of approximately 32 kDa (predicted molecular weight: 32 kDa). Unpurified format: 1/5,000 - 1/20,000
IHC-P		1/50 - 1/500. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.
IP		1/50.

ターゲット情報

機能	Inhibitor of protein-phosphatase 1.
配列類似性	Belongs to the protein phosphatase inhibitor 1 family.
翻訳後修飾	Dopamine- and cyclic AMP-regulated neuronal phosphoprotein. Phosphorylation of Thr-34 is required for activity.
細胞内局在	Cytoplasm.

画像



Western blot - Anti-DARPP32 antibody [EP721Y] (ab40802)


Anti-DARPP32 antibody [EP721Y] (ab40802) at 1/50000 dilution + Human brain tissue lysates. at 10 µg

Secondary

Anti-Rabbit IgG (HRP), specific to the non-reduced form of IgG at 1/10000 dilution

Predicted band size: 32 kDa


Additional bands at: 32 kDa. We are unsure as to the identity of these extra bands.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-DARPP32 antibody [EP721Y] (ab40802)

ab40802 staining DARPP32 in human cerebral cortex tissue sections by Immunohistochemistry (IHC-P - paraformaldehyde-fixed, paraffin-embedded sections). Antigen retrieval was by heat mediation in a EDTA buffer. Samples were incubated with primary antibody at a dilution of 1/500. A goat anti-rabbit IgG H&L (HRP) **ab97051** was used as the secondary antibody.

Negative control 1: PBS in place of primary antibody.




Immunocytochemistry/ Immunofluorescence - Anti-DARPP32 antibody [EP721Y] (ab40802)

ab40802 staining DARPP32 in SH-SY5Y (human neuroblastoma) cells by ICC/IF (Immunocytochemistry/immunofluorescence). Cells were fixed with 100% methanol and permeabilized with 0.1% Triton X-100. Samples were incubated with primary antibody at a dilution of 1/100. A goat anti rabbit IgG (Alexa Fluor® 488) (**ab150077**) was used as the secondary antibody at a dilution of 1/1000. **ab7291** and **ab150120** were used as counterstains for primary antibody ab40802 and secondary antibody **ab150077** respectively and DAPI was used as a nuclear counterstain.


Negative control 1: Rabbit primary antibody and anti-mouse secondary antibody (**ab150120**)

Negative control 2: Mouse primary antibody (**ab7291**) and anti-rabbit secondary antibody (**ab150077**)



Immunohistochemical analysis human gastric adenocarcinoma using ab40802. (dilution 1:250)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-DARPP32 antibody [EP721Y] (ab40802)



ab40802 immunoprecipitating DARPP32. 10µg of cell lysate was incubated with primary antibody at a dilution of 1/50 and VeriBlot for IP Detection Reagent (HRP) ([ab131366](#)) at a dilution of 1/10000.


Lane 1: Human fetal forebrain whole cell lysate (10ug)

Lane 2: Human fetal forebrain whole cell lysate

Lane 3: Rabbit monoclonal IgG ([ab172730](#)) instead of ab40802 in Human fetal forebrain whole cell lysate

Immunoprecipitation - Anti-DARPP32 antibody [EP721Y] (ab40802)

Why choose a recombinant antibody?



Anti-DARPP32 antibody [EP721Y] (ab40802)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors