

Anti-DFNA5/GSDME antibody [EPR19859-60] - C-terminal ab221843

KO 評価済 リコンビナント RabMAb

1 References 画像数 5

製品の概要

製品名	Anti-DFNA5/GSDME antibody [EPR19859-60] - C-terminal
製品の詳細	Rabbit monoclonal [EPR19859-60] to DFNA5/GSDME - C-terminal
由来種	Rabbit
特異性	Expression levels of the target protein vary with sample type and some optimisation may be required.
アプリケーション	適用あり: WB, IP
種交差性	交差種: Human
免疫原	Recombinant full length protein. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: SH-SY5Y (human neuroblastoma epithelial cell), whole cell lysate. HEK-293 transfected with a myc tagged, DFNA5/GSDME (C-terminus) expression vector, whole cell lysate. Human fetal brain lysate. IP: SH-SY5Y (human neuroblastoma epithelial cell), whole cell lysate.
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none">- High batch-to-batch consistency and reproducibility- Improved sensitivity and specificity- Long-term security of supply- Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Preservative: 0.01% Sodium azide Constituents: 0.05% BSA, 40% Glycerol (glycerin, glycerine), PBS

精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR19859-60
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab221843の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		1/1000. Detects a band of approximately 36, 55 kDa (predicted molecular weight: 54 kDa).
IP		1/40.

ターゲット情報

組織特異性

Expressed in cochlea. Low level of expression in heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas, with highest expression in placenta.


関連疾患

Defects in DFNA5 are the cause of deafness autosomal dominant type 5 (DFNA5) [MIM:600994]. DFNA5 is a form of sensorineural hearing loss. Sensorineural deafness results from damage to the neural receptors of the inner ear, the nerve pathways to the brain, or the area of the brain that receives sound information.

配列類似性

Belongs to the gasdermin family.

画像



All lanes : Anti-DFNA5/GSDME antibody [EPR19859-60] - C-terminal (ab221843) at 1/5000 dilution

Lane 1 : SH-SY5Y (human neuroblastoma epithelial cell), whole cell lysate

Lane 2 : DFNA5/GSDME knockout SH-SY5Y, whole cell lysate

Lysates/proteins at 10 µg per lane.

Secondary


All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/100000 dilution

Predicted band size: 54 kDa

Exposure time: 30 seconds

Blocking/Dilution: 5% NFDm/TBST.

The cell lysates were kindly provided by our collaborator, Dr. Feng Shao, NIBS.



Western blot - Anti-DFNA5/GSDME antibody [EPR19859-60] - C-terminal (ab221843)

All lanes : Anti-DFNA5/GSDME antibody [EPR19859-60] - C-terminal (ab221843) at 1/1000 dilution

Lane 1 : HEK-293 (human embryonic kidney epithelial cell) transfected with a control vector containing DDDDK tag, whole cell lysate

Lane 2 : HEK-293 transfected with a DDDDK tagged, DFNA5/GSDME (N-terminus) expression vector, whole cell lysate

Lane 3 : HEK-293 transfected with a myc tagged, DFNA5/GSDME (C-terminus) expression vector, whole cell lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/100000 dilution

Developed using the ECL technique.


Predicted band size: 54 kDa

Exposure time: 5 seconds

Blocking/Dilution: 5% NFDm/TBST .

This antibody is specific to the C-terminus of DFNA5/GSDME.

The cell lysates were kindly provided by our collaborator, Dr. Feng Shao, NIBS.



Western blot - Anti-DFNA5/GSDME antibody
[EPR19859-60] - C-terminal (ab221843)

Anti-DFNA5/GSDME antibody [EPR19859-60] - C-terminal
(ab221843) at 1/1000 dilution + human fetal brain at 10 µg

Secondary

VeriBlot for IP Detection Reagent (HRP) ([ab131366](#)) at 1/4000
dilution


Developed using the ECL technique.

Predicted band size: 54 kDa

Exposure time: 3 minutes

Blocking/Dilution: 5% NFDm/TBST.

The molecular weight observed is consistent with what has been
described in the literature (PMID: 28459430).



Immunoprecipitation - Anti-DFNA5/GSDME antibody
[EPR19859-60] - C-terminal (ab221843)

DFNA5/GSDME was immunoprecipitated from 0.35 mg of SH-
SY5Y (human neuroblastoma epithelial cell) whole cell lysate with
ab221843 at 1/40 dilution. Western blot was performed from the
immunoprecipitate using ab221843 at 1/500 dilution. VeriBlot for IP
Detection Reagent (HRP) ([ab131366](#)), was used for detection at
1/1000 dilution.

Lane 1: SH-SY5Y whole cell lysate 10ug (Input).

Lane 2: ab221843 IP in SH-SY5Y whole cell lysate.

Lane 3: Rabbit monoclonal IgG ([ab172730](#)) instead of ab221843 in
SH-SY5Y whole cell lysate.

Blocking and dilution buffer and concentration: 5% NFDm/TBST.

Observed MW 55 kDa.

Exposure time: 1 second.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-DFNA5/GSDME antibody [EPR19859-60] - C-terminal (ab221843)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors