

Anti-DHX58/RLR antibody ab67272

★★★★★ [1 Abreviews](#) [画像数 1](#)

製品の概要

製品名	Anti-DHX58/RLR antibody
製品の詳細	Rabbit polyclonal to DHX58/RLR
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Rat
免疫原	Synthetic peptide corresponding to Human DHX58/RLR (internal sequence).
ポジティブ・コントロール	WB: Rat kidney tissue lysate.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab67272の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (1)	Use a concentration of 1 - 2 µg/ml. Detects a band of approximately 45, 77 kDa (predicted molecular weight: 45, 77 kDa).

ターゲット情報

機能	Participates in innate immune defense against viruses. Upon interaction with intracellular dsRNA produced during viral replication, triggers a transduction cascade which results in the induction and expression of antiviral cytokines such as IFN-beta and RANTES (CCL5). The RNA helicase domain may recognize and structurally modify viral RNA to facilitate detection by DDX58/RIG-I or by IFIH1/MDA5 whose affinity for dsRNA is lower.
配列類似性	Belongs to the helicase family. Contains 1 helicase ATP-binding domain. Contains 1 helicase C-terminal domain.
細胞内局在	Cytoplasm.

画像



Lane 1 : Anti-DHX58/RLR antibody (ab67272) at 1 µg/ml

Lane 2 : Anti-DHX58/RLR antibody (ab67272) at 2 µg/ml

All lanes : Rat kidney tissue lysate

Lysates/proteins at 15 µg per lane.

Secondary

All lanes : anti rabbit secondary antibody

Predicted band size: 45, 77 kDa

Observed band size: 45,77 kDa

There are 4 isoforms of LGP2 - 376, 678, 379, and 368 amino acids. The larger band is isoform b and the smaller one is either a, c, or d since they are similar in size.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors