

Anti-DR5 antibody [EPR19310] ab199357

KO 評価済 リコンビナント RabMAb

★★★★★ [5 Abreviews](#) [12 References](#) [画像数 10](#)

製品の概要

製品名	Anti-DR5 antibody [EPR19310]
製品の詳細	Rabbit monoclonal [EPR19310] to DR5
由来種	Rabbit
アプリケーション	適用あり: Flow Cyt (Intra), WB, ICC/IF, IP
種交差性	交差種: Human
免疫原	Recombinant fragment. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: Untreated and treated with 0.5ÅµM/ml doxorubicin for 24 hours HCT 116 whole cell lysates; HeLa, HAP1, HepG2 and HT1080 whole cell lysates; Human melanoma lysate. ICC/IF: HT1080 and HCT 116 cells. Flow Cyt (intra): HCT 116 cells. IP: HT-1080 treated with 5ÅµM MG132 for 4 hour whole cell lysate; HCT 116 whole cell lysate.
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb® technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb® patents.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.2</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 59% PBS, 40% Glycerol, 0.05% BSA</p>
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR19310

アイソタイプ

IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab199357の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt (Intra)		1/70.
WB	★★★★★ (2)	1/1000. Detects a band of approximately 48, 40 kDa (predicted molecular weight: 48 kDa).
ICC/IF		1/100.
IP		1/30.

ターゲット情報

機能	Receptor for the cytotoxic ligand TNFSF10/TRAIL. The adapter molecule FADD recruits caspase-8 to the activated receptor. The resulting death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation which initiates the subsequent cascade of caspases (aspartate-specific cysteine proteases) mediating apoptosis. Promotes the activation of NF-kappa-B. Essential for ER stress-induced apoptosis.
組織特異性	Widely expressed in adult and fetal tissues; very highly expressed in tumor cell lines such as HeLaS3, K-562, HL-60, SW480, A-549 and G-361; highly expressed in heart, peripheral blood lymphocytes, liver, pancreas, spleen, thymus, prostate, ovary, uterus, placenta, testis, esophagus, stomach and throughout the intestinal tract; not detectable in brain.
関連疾患	Squamous cell carcinoma of the head and neck
配列類似性	Contains 1 death domain. Contains 3 TNFR-Cys repeats.
細胞内局在	Membrane.

画像



Western blot - Anti-DR5 antibody [EPR19310]
(ab199357)

All lanes : Anti-DR5 antibody [EPR19310] (ab199357) at 1/1000 dilution

Lane 1 : Wild-type HeLa cell lysate

Lane 2 : DR5 knockout HeLa cell lysate

Lysates/proteins at 20 µg per lane.


Performed under reducing conditions.

Predicted band size: 48 kDa

Observed band size: 47 kDa


Lanes 1- 2: Merged signal (red and green). Green - ab199357 observed at 47 kDa. Red - Anti-GAPDH antibody [6C5] - Loading Control ([ab8245](#)) observed at 37 kDa.

ab199357 was shown to react with DR5 in wild-type HeLa cells in western blot. Loss of signal was observed when knockout cell line [ab264922](#) (knockout cell lysate [ab257748](#)) was used. Wild-type HeLa and DR5 knockout HeLa cell lysates were subjected to SDS-PAGE. Membrane was blocked for 1 hour at room temperature in 0.1% TBST with 3% non-fat dried milk. ab199357 and Anti-GAPDH antibody [6C5] - Loading Control ([ab8245](#)) overnight at 4°C at a 1 in 1000 dilution and a 1 in 20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye®800CW) preadsorbed ([ab216773](#)) and Goat anti-Mouse IgG H&L (IRDye®680RD) preadsorbed ([ab216776](#)) secondary antibodies at 1 in 20000 dilution for 1 hour at room temperature before imaging.



Flow Cytometry (Intracellular) - Anti-DR5 antibody [EPR19310] (ab199357)

Intracellular flow cytometric analysis of 4% paraformaldehyde-fixed HCT 116 (Human colorectal carcinoma cell line) cells labeling DR5 with ab199357 at 1/70 dilution (red) compared with a Rabbit IgG, monoclonal [EPR25A] - Isotype Control (**ab172730**) (black) and an unlabelled control (cells without incubation with primary antibody and secondary antibody; blue). Goat anti Rabbit IgG (Alexa Fluor® 488) at 1/2000 dilution was used as the secondary antibody.



Immunocytochemistry/ Immunofluorescence - Anti-DR5 antibody [EPR19310] (ab199357)


Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized HT1080 (Human fibrosarcoma cell line) cells labeling DR5 with ab199357 at 1/100 dilution, followed by Goat Anti-Rabbit IgG (Alexa Fluor® 488) (**ab150077**) secondary antibody at 1/1000 dilution (green). Confocal image showing cytoplasmic staining on HT1080 cells. The nuclear counter stain is DAPI (blue).

Tubulin is detected with Anti-alpha Tubulin antibody [DM1A] - Loading Control (**ab7291**) at 1/1000 dilution and Goat Anti-Mouse IgG H&L (Alexa Fluor® 594) preadsorbed (**ab150120**) at 1/1000 dilution (red).

The negative controls are as follows:

-ve control 1: ab199357 at 1/100 dilution followed by **ab150120** at 1/1000 dilution.

-ve control 2: **ab7291** at 1/1000 dilution followed by **ab150077** at 1/1000 dilution.



Western blot - Anti-DR5 antibody [EPR19310]
(ab199357)

All lanes : Anti-DR5 antibody [EPR19310] (ab199357) at 1/1000 dilution

Lane 1 : Wild-type HAP1 whole cell lysate

Lane 2 : TNFRSF10B knockout HAP1 whole cell lysate


Lysates/proteins at 20 µg per lane.

Predicted band size: 48 kDa

Observed band size: 47 kDa

Lanes 1 - 2: Merged signal (red and green). Green - ab199357 observed at 47 kDa. Red - loading control, **ab9484**, observed at 37 kDa.

ab199357 was shown to recognize DR5 in wild-type HAP1 cells as signal was lost at the expected MW in TNFRSF10B knockout cells. Additional cross-reactive bands were observed in the wild-type and knockout cells. Wild-type and TNFRSF10B knockout samples were subjected to SDS-PAGE. Ab199357 and **ab9484** (Mouse anti-GAPDH loading control) were incubated overnight at 4°C at 1/1000 dilution and 1/20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preabsorbed **ab216773** and Goat anti-Mouse IgG H&L (IRDye® 680RD) preabsorbed **ab216776** secondary antibodies at 1/20000 dilution for 1 hour at room temperature before imaging.



Western blot - Anti-DR5 antibody [EPR19310]
(ab199357)

All lanes : Anti-DR5 antibody [EPR19310] (ab199357) at 1/1000 dilution

Lane 1 : Untreated HCT 116 (Human colorectal carcinoma cell line) whole cell lysate

Lane 2 : HCT 116 (Human colorectal carcinoma cell line) treated with 0.5 μM/ml doxorubicin for 24 hours whole cell lysate

Lysates/proteins at 20 μg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/100000 dilution


Predicted band size: 48 kDa

Observed band size: 40,48 kDa

Exposure time: 8 seconds

Blocking/Dilution buffer: 5% NFDm/TBST.

Doxorubicin treatment elevated the expression of DR5 (PMID: 12496481; PMID: 11468181; PMID: 11090076). The expression profile is consistent with what has been described in the literature (PMID:20515924; PMID:16297203).



All lanes : Anti-DR5 antibody [EPR19310] (ab199357) at 1/1000 dilution

Lane 1 : Human melanoma lysate

Lane 2 : HeLa (Human epithelial cell line from cervix adenocarcinoma) whole cell lysate

Lane 3 : HepG2 (Human liver hepatocellular carcinoma cell line) whole cell lysate

Lane 4 : HT1080 (Human fibrosarcoma cell line) whole cell lysate

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/100000 dilution


Predicted band size: 48 kDa

Observed band size: 40,48 kDa

Blocking/Dilution buffer: 5% NFDM/TBST.

Exposure time: Lane 1-2: 3 minutes; Lane 3: 10 seconds; Lane 4: 8 seconds.

The expression profile is consistent with what has been described in the literature (PMID:20515924; PMID:16297203).




Immunocytochemistry/ Immunofluorescence - Anti-DR5 antibody [EPR19310] (ab199357)

Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized HCT 116 (Human colorectal carcinoma cell line) cells labeling DR5 with ab199357 at 1/100 dilution, followed by Goat Anti-Rabbit IgG (Alexa Fluor® 488) (**ab150077**) secondary antibody at 1/1000 dilution (green). Confocal image showing membranous and cytoplasmic staining on HCT 116 cells. The nuclear counter stain is DAPI (blue). Tubulin is detected with Anti-alpha Tubulin antibody [DM1A] - Loading Control (**ab7291**) at 1/1000 dilution and Goat Anti-Mouse IgG H&L (Alexa Fluor® 594) preadsorbed (**ab150120**) at 1/1000 dilution (red).

The negative controls are as follows:

-ve control 1: ab199357 at 1/100 dilution followed by **ab150120** at 1/1000 dilution.

-ve control 2: **ab7291** at 1/1000 dilution followed by **ab150077** at 1/1000 dilution.



Immunoprecipitation - Anti-DR5 antibody [EPR19310] (ab199357)

DR5 was immunoprecipitated from 0.35mg of HT1080 (Human fibrosarcoma cell line) treated with 5μM MG132 for 4 hour whole cell lysate with ab199357 at 1/30 dilution. Western blot was performed from the immunoprecipitate using ab199357 at 1/1000 dilution. VeriBlot for IP Detection Reagent (HRP) (**ab131366**), was used for detection at 1/10000 dilution.


Lane 1: HT1080 treated with 5μM MG132 for 4 hour whole cell lysate, 10μg (Input).

Lane 2: ab199357 IP in HT1080 treated with 5μM MG132 for 4 hour whole cell lysate.

Lane 3: Rabbit IgG, monoclonal [EPR25A] - Isotype Control (**ab172730**) instead of ab199357 in HT1080 treated with 5μM MG132 for 4 hour whole cell lysate.

Blocking and dilution buffer and concentration: 5% NFDm/TBST.

Exposure time: 10 seconds.



Immunoprecipitation - Anti-DR5 antibody
[EPR19310] (ab199357)

DR5 was immunoprecipitated from 0.35mg of HCT 116 (Human colorectal carcinoma cell line) whole cell lysate with ab199357 at 1/30 dilution. Western blot was performed from the immunoprecipitate using ab199357 at 1/1000 dilution. VeriBlot for IP Detection Reagent (HRP) ([ab131366](#)), was used for detection at 1/10000 dilution.

Lane 1: HCT 116 whole cell lysate, 10µg (Input).

Lane 2: ab199357 IP in HCT 116 whole cell lysate.

Lane 3: Rabbit IgG, monoclonal [EPR25A] - Isotype Control ([ab172730](#)) instead of ab199357 in HCT 116 whole cell lysate.

Blocking and dilution buffer and concentration: 5% NFDm/TBST.

Exposure time: 10 seconds.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-DR5 antibody [EPR19310] (ab199357)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors