

FITC Anti-C Reactive Protein antibody ab34659

3 References

医薬用外劇物

製品の概要

製品名	FITC Anti-C Reactive Protein antibody
製品の詳細	FITC Goat polyclonal to C Reactive Protein
由来種	Goat
標識	FITC. Ex: 493nm, Em: 528nm
アプリケーション	適用あり: Immunomicroscopy, Flow Cyt
種交差性	交差種: Human
免疫原	Full length protein corresponding to Human C Reactive Protein. Database link: P02741
ポジティブ・コントロール	Partially purified human CRP
特記事項	<p>2.2moles Fluorescein (Diaminotriazynlamino-fluorescein (DTAF0), 530 daltons) per mole of Goat IgG. Absorption Wavelength: 495 nm Emission Wavelength: 528 nm.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Do Not Freeze.
バッファー	<p>pH: 6.50</p> <p>Preservative: 0.01% Thimerosal (merthiolate)</p> <p>Constituents: 1% BSA, 0.812% Sodium chloride, 0.164% Sodium phosphate</p> <p>BSA is IgG and Protease free</p>
精製度	IgG fraction
特記事項 (精製)	This product is an IgG fraction antibody purified from monospecific antiserum by a multi-step

process including delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against buffer.

ポリ/モノ

ポリクローナル

アイソタイプ

IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab34659の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Immunomicroscopy		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. ab37394 - Goat polyclonal IgG, is suitable for use as an isotype control with this antibody.

ターゲット情報

機能	Displays several functions associated with host defense: it promotes agglutination, bacterial capsular swelling, phagocytosis and complement fixation through its calcium-dependent binding to phosphorylcholine. Can interact with DNA and histones and may scavenge nuclear material released from damaged circulating cells.
組織特異性	Found in plasma.
配列類似性	Belongs to the pentaxin family. Contains 1 pentaxin domain.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors