

Anti-Fumonisin antibody [2A2] ab42124

1 References

製品の概要

製品名	Anti-Fumonisin antibody [2A2]
製品の詳細	Mouse monoclonal [2A2] to Fumonisin
由来種	Mouse
アプリケーション	適用あり: ELISA
種交差性	交差種: Fusarium verticillioides
免疫原	Chemical/ Small Molecule conjugated to BSA: Fumonisin
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Constituents: 0.75% Glycine, 1.21% Tris, 2% Sucrose
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	2A2
アイソタイプ	IgG1
軽鎖の種類	lambda

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab42124の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ELISA		1/2000 - 1/10000.

ターゲット情報

関連性 Fumonisin is a group of toxins (B1 or FB1, B2 or FB2, which are equally toxic, and B3 or FB3, which is relatively nontoxic) which are produced by fungi (*Fusarium moniliforme* and *F. proliferatum*) that causes equine leukoencephalomalacia, a disease of the central nervous system of equines and is associated with the feeding of moldy corn. Swine may develop porcine pulmonary edema.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors