

Anti-GAL4 antibody [EPR12710(B)] ab175185

リコンビナント **RabMAb**

画像数 4

製品の概要

製品名	Anti-GAL4 antibody [EPR12710(B)]
製品の詳細	Rabbit monoclonal [EPR12710(B)] to GAL4
由来種	Rabbit
アプリケーション	適用あり: IHC-P, WB 適用なし: Flow Cyt, ICC/IF or IP
種交差性	交差種: Rat, Human 非交差種: Mouse
免疫原	Synthetic peptide within Human GAL4 aa 1-100 (Cysteine residue). The exact sequence is proprietary. Database link: P56470
ポジティブ・コントロール	HT-29, Human small intestine and rat small intestine lysates. Human colon and colon carcinoma tissue.
特記事項	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information see here . Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 9% PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA, 50% Tissue culture supernatant


アプリケーション

アプリケーション	Abreviews	特記事項
IHC-P		1/250 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
WB		1/1000 - 1/5000. Predicted molecular weight: 36 kDa.


ターゲット情報

ドメイン Contains two homologous but distinct carbohydrate-binding domains.

画像



Predicted band size: 36 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-GAL4 antibody
[EPR12710(B)] (ab175185)

Immunohistochemical analysis of paraffin-embedded Human colon carcinoma tissue labeling GAL4 with ab175185 at a 1/250 dilution.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-GAL4 antibody
[EPR12710(B)] (ab175185)

Immunohistochemical analysis of paraffin-embedded Human colon tissue labeling GAL4 with ab175185 at a 1/250 dilution.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-GAL4 antibody [EPR12710(B)] (ab175185)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors