


Anti-Gata6 antibody ab22600

★★★★☆ [5 Abreviews](#) [28 References](#) [画像数 1](#)

製品の概要

製品名	Anti-Gata6 antibody
製品の詳細	Rabbit polyclonal to Gata6
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse 交差が予測される動物種: Rat, Chicken, Cow, Human, Pig 
免疫原	Synthetic peptide corresponding to Human Gata6 aa 400 to the C-terminus conjugated to keyhole limpet haemocyanin. (Peptide available as ab24410 , ab24411)
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Preservative: 0.02% Sodium azide Constituent: PBS Batches of this product that have a concentration < 1mg/ml may have BSA added as a stabilising agent. If you would like information about the formulation of a specific lot, please contact our scientific support team who will be happy to help.
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab22600の使用に適用されます**


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★★ (3)	Use a concentration of 1 µg/ml. Detects a band of approximately 45 kDa (predicted molecular weight: 45 kDa). Abcam recommends using milk as the blocking agent.

ターゲット情報

機能	Thought to be important for regulating terminal differentiation and/or proliferation.
組織特異性	Expressed in heart, gut and gut-derived tissues.
配列類似性	Contains 2 GATA-type zinc fingers.
ドメイン	The GATA-type zinc fingers mediate interaction with LMCD1.
細胞内局在	Nucleus.

画像



Anti-Gata6 antibody (ab22600) at 1 µg/ml + D3 ES Cells Day 15
EB (Gata6 +ve) at 10 µg

Secondary

IRDye 680 Conjugated Goat Anti-Rabbit IgG at 1/10000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Predicted band size: 45 kDa

Observed band size: 45 kDa

Additional bands at: 73 kDa. We are unsure as to the identity of these extra bands.

Exposure time: 2 minutes

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors