

Anti-GFAP antibody [RM1003] ab278054

リコンビナント RabMAb

画像数 12

製品の概要

製品名	Anti-GFAP antibody [RM1003]
製品の詳細	Rabbit recombinant multiclonal [RM1003] to GFAP
由来種	Rabbit
アプリケーション	適用あり: Flow Cyt (Intra), WB, IHC-P, IHC-Fr, IP, ICC/IF
種交差性	交差種: Mouse, Rat, Human
免疫原	Recombinant fragment. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: Human, mouse and rat brain tissue lysates; U-87 MG, C6 whole cell lysates. IHC-P: Human, mouse and rat cerebrum tissue. IHC-Fr: Mouse cerebrum tissue. Flow Cyt: Mouse primary brain cells. IP: Mouse brain lysate. ICC/IF: Rat hippocampal mixed glia and mouse primary neural/glia cells.
特記事項	Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.01% Sodium azide Constituents: 59.94% PBS, 0.05% BSA, 40% Glycerol (glycerin, glycerine)
精製度	Protein A purified
ポリ/モノ	Recombinant Multiclonal
クローン名	RM1003
アイソタイプ	IgG


アプリケーション

The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab278054の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Flow Cyt (Intra)		1/500.
WB		1/1000. Detects a band of approximately 40-54 kDa (predicted molecular weight: 49 kDa).
IHC-P		1/2000. Perform heat mediated antigen retrieval with Tris/EDTA buffer pH 9.0 before commencing with IHC staining protocol.
IHC-Fr		1/700. Heat mediated antigen retrieval using sodium citrate buffer (10mM citrate pH 6.0 + 0.05% Tween-20).
IP		1/30.
ICC/IF		Use a concentration of 0.1 - 1 µg/ml.

ターゲット情報	
機能	GFAP, a class-III intermediate filament, is a cell-specific marker that, during the development of the central nervous system, distinguishes astrocytes from other glial cells.
組織特異性	Expressed in cells lacking fibronectin.
関連疾患	Defects in GFAP are a cause of Alexander disease (ALEXD) [MIM:203450]. Alexander disease is a rare disorder of the central nervous system. It is a progressive leukoencephalopathy whose hallmark is the widespread accumulation of Rosenthal fibers which are cytoplasmic inclusions in astrocytes. The most common form affects infants and young children, and is characterized by progressive failure of central myelination, usually leading to death usually within the first decade. Infants with Alexander disease develop a leukoencephalopathy with macrocephaly, seizures, and psychomotor retardation. Patients with juvenile or adult forms typically experience ataxia, bulbar signs and spasticity, and a more slowly progressive course.
配列類似性	Belongs to the intermediate filament family.
翻訳後修飾	Phosphorylated by PKN1.
細胞内局在	Cytoplasm. Associated with intermediate filaments.

画像



Immunocytochemistry/ Immunofluorescence - Anti-GFAP antibody [RM1003] (ab278054)


Immunofluorescent analysis of 4% Paraformaldehyde-fixed, 0.1% TritonX-100 permeabilized mouse primary neural/glia cells labelling GFAP with ab278054 at 1/500 (0.938 µg/ml) dilution, followed by **ab150077** Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) antibody at 1/1000 dilution (Green).

Confocal image showing cytoplasmic staining in mouse primary astrocytes. Confocal scanning Z step was set as 0.3 µm followed by image processing with maximum Z projection.

ab10062 Anti-GFAP mouse monoclonal antibody at 1/200 dilution, followed by **ab150120** Goat Anti-Mouse IgG H&L (Alexa Fluor® 594) at a 1/500 dilution (Red) to counterstain. The nuclear counterstain was DAPI (Blue).

Negative control 1: ab278054 at a 1/500 dilution followed by **ab150120** at a 1/500 dilution.

Negative control 2: **ab10062** at a 1/500 dilution followed by **ab150077** at a 1/1000 dilution.




Immunocytochemistry/ Immunofluorescence - Anti-GFAP antibody [RM1003] (ab278054)

Immunofluorescence staining of GFAP using ab278054 in primary rat hippocampal mixed glia, (prepared from P2 rat hippocampal brain area, obtained from Transnetyx Tissue by BrainBits, LLC, cat.no. SDPHP4m), DIV4. The cells were fixed with 100% MeOH (5 min), permeabilized with 0.1% Triton-X-100 (in PBS) for 5 mins and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Tween for 1h. The cells were then incubated overnight at +4°C with ab278054 at 0.1 µg/ml and **ab4674**, Anti-GFAP antibody, at 1/1000 dilution. Cells were then incubated with **ab150081**, Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) preadsorbed at 1/1000 dilution (shown in green) and **ab150176**, Goat Anti-Chicken IgY H&L (Alexa Fluor® 594) preadsorbed at 1/1000 dilution (shown in red). Nuclear DNA was labelled with DAPI (shown in blue).

Images were acquired with the Perkin Elmer Operetta HCA and a maximum intensity projection of confocal sections is shown. The antibody ab278054 gave comparable results using 4%

formaldehyde fixation (10 min).



All lanes : Anti-GFAP antibody [RM1003] (ab278054) at 1/1000 dilution

Lane 1 : U-87 MG (Human glioblastoma-astrocytoma epithelial cell) whole cell lysate

Lane 2 : SH-SY5Y (Human neuroblastoma epithelial cell) whole cell lysate

Lane 3 : SK-N-BE (2) (Human neuroblastoma neuroblast) whole cell lysate

Lane 4 : IMR-32 (Human lung fibroblast) whole cell lysate

Lane 5 : HeLa (Human cervix adenocarcinoma epithelial cell) whole cell lysate

Lane 6 : Neuro-2a (Mouse neuroblastoma neuroblast) whole cell lysate

Lane 7 : NIH/3T3 (Mouse embryonic fibroblast) whole cell lysate

Lane 8 : C6 (Rat glial tumor glial cell) whole cell lysate

Lane 9 : PC-12 (Rat adrenal gland pheochromocytoma) whole cell lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG (HRP), specific to the non-reduced form of IgG at 1/20000 dilution

Predicted band size: 49 kDa

Observed band size: 40-54 kDa

Exposure time: 180 seconds


Blocking and Diluting buffer and concentration: 5% NFDM/TBST

GFAP is expressed primarily by astrocytes in the CNS as astrocytes marker. We have no suitable astrocyte for western blot testing.

U-87 MG and C6 were reported to express low level of GFAP (PMID: 23839947, PMID: 22096544, PMID: 20669222).

NIH/3T3, IMR-32, HeLa and PC-12 were reported as negative cell lines for GFAP (PMID: 824020; PMID:2294; PMID: 6340792; PMID:9466565, PMID:7895062, PMID: 28700643, PMID: 19272755).

No literature was found to support the expression in SH-SY5Y, Neuro-2a and SK-N-BE (2) cell.




Immunohistochemistry (Frozen sections) - Anti-GFAP antibody [RP1003] (ab278054)

Immunohistochemical analysis of 4% PFA-fixed, 0.2% Triton X-100 permeabilized mouse cerebrum tissue labeling GFAP with ab278054 at 1/700 (0.938 µg/ml) dilution followed by **ab150077** Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) at a 1/1000 dilution. Nuclear counterstain is DAPI. Positive staining on mouse cerebrum.

Secondary antibody only control: Secondary antibody is **ab150077** Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) used at a 1/1000 dilution.

Heat mediated antigen retrieval using sodium citrate buffer (10 mM citrate pH 6.0 + 0.05% Tween-20).




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-GFAP antibody [RP1003] (ab278054)


Immunohistochemical analysis of paraffin-embedded human cerebrum tissue labeling GFAP with ab278054 at 1/2000 (0.235 µg/ml) dilution followed by ready to use Rabbit specific IHC polymer detection kit HRP/DAB (**ab209101**). Positive staining on human cerebrum. The section was incubated with ab278054 for 30 mins at room temperature. The immunostaining was performed on a Leica Biosystems BOND® RX instrument. Counterstained with Hematoxylin.

Secondary antibody only control: Secondary antibody is a ready to use Rabbit specific IHC polymer detection kit HRP/DAB (**ab209101**).

Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0, epitope retrieval solution2) for 20 mins.



Flow Cytometry (Intracellular) - Anti-GFAP antibody [RP1003] (ab278054)



All lanes : Anti-GFAP antibody [RM1003] (ab278054) at 1/1000 dilution

Lane 1 : Mouse brain lysate

Lane 2 : Rat brain lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/100000 dilution

Predicted band size: 49 kDa


Observed band size: 40-54 kDa

Exposure time: 3 seconds

The molecular weight observed is consistent with the literature

(PMID: 25975510).

Blocking/Dilution buffer: 5% NFDM/TBST.



GFAP was immunoprecipitated from 0.35 mg mouse brain lysate with ab278054 at 1/30 dilution (2µg in 0.35mg lysates). Western blot was performed on the immunoprecipitate using ab278054 at 1/1000 dilution. VeriBlot for IP Detection Reagent (HRP) ([ab131366](#)) was used at 1/5000 dilution.

Lane 1: Mouse brain lysate 10 µg.


Lane 2: ab278054 IP in mouse brain lysate.

Lane 3: Rabbit monoclonal IgG ([ab172730](#)) instead of ab278054 in mouse brain lysate.

Blocking and dilution buffer and concentration: 5% NFDM/TBST.

Exposure time: 3 seconds.


The molecular weight observed is consistent with the literature (PMID: 824020; PMID:2294; PMID: 6340792).



Immunohistochemical analysis of paraffin-embedded mouse cerebrum tissue labeling GFAP with ab278054 at 1/2000 (0.235 µg/ml) dilution followed by ready to use Rabbit specific IHC polymer detection kit HRP/DAB ([ab209101](#)). Positive staining on mouse cerebrum. The section was incubated with ab278054 for 30 mins at room temperature. The immunostaining was performed on a Leica Biosystems BOND® RX instrument. Counterstained with Hematoxylin.

Secondary antibody only control: Secondary antibody is a ready to use Rabbit specific IHC polymer detection kit HRP/DAB ([ab209101](#)).

Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0, epitope retrieval solution2) for 20 mins.



Western blot - Anti-GFAP antibody [RP1003]
(ab278054)

All lanes : Anti-GFAP antibody [RM1003] (ab278054) at 1/1000 dilution

Lane 1 : Human brain lysate

Lane 2 : Human liver lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/100000 dilution

Predicted band size: 49 kDa


Observed band size: 40-54 kDa

Exposure time: 3 seconds

Negative control: Human liver. (PMID: 25975510).

The molecular weight observed is consistent with the literature (PMID: 25975510).

Blocking/Diluting buffer: 5% NFDm/TBST.







Immunohistochemical analysis of paraffin-embedded rat cerebrum tissue labeling GFAP with ab278054 at 1/2000 (0.235 µg/ml) dilution followed by ready to use Rabbit specific IHC polymer detection kit HRP/DAB (**ab209101**). Positive staining on rat cerebrum. The section was incubated with ab278054 for 30 mins at room temperature. The immunostaining was performed on a Leica Biosystems BOND® RX instrument. Counterstained with Hematoxylin.

Secondary antibody only control: Secondary antibody is a ready to use Rabbit specific IHC polymer detection kit HRP/DAB (**ab209101**).

Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0, epitope retrieval solution2) for 20 mins.

Why choose a recombinant antibody?

 <p>Research with confidence Consistent and reproducible results</p>	 <p>Long-term and scalable supply Recombinant technology</p>
 <p>Success from the first experiment Confirmed specificity</p>	 <p>Ethical standards compliant Animal-free production</p>

Anti-GFAP antibody [RP1003] (ab278054)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors