

Anti-GM-CSF antibody ab9667

[2 References](#) [画像数 1](#)

製品の概要

製品名	Anti-GM-CSF antibody
製品の詳細	Rabbit polyclonal to GM-CSF
由来種	Rabbit
アプリケーション	適用あり: Sandwich ELISA
種交差性	交差種: Recombinant fragment
免疫原	Recombinant full length protein corresponding to Human GM-CSF. Highly pure (>98%) recombinant hGM-CSF (human Granulocyte Macrophage Colony Stimulating Factor) Database link: P04141
ポジティブ・コントロール	Recombinant human GM-CSF protein (ab9668) can be used as a positive control in WB.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Lyophilized:Reconstitute with 200µl of sterile water. Please note that if you receive this product in liquid form it has already been reconstituted as described and no further reconstitution is necessary.
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C.
バッファー	No preservative, sterile filtered
キャリア・フリー	はい
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab9667の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Sandwich ELISA		Use a concentration of 0.5 - 2 µg/ml. Allows the detection of 0.2 - 0.4 ng/well of recombinant hGM-CSF (using 100µl/well antibody solution). Can be paired for Sandwich ELISA with ab84259 . Can be used as the capture antibody when paired with recommended pair.

ターゲット情報

機能	Cytokine that stimulates the growth and differentiation of hematopoietic precursor cells from various lineages, including granulocytes, macrophages, eosinophils and erythrocytes.
配列類似性	Belongs to the GM-CSF family.
細胞内局在	Secreted.

画像



Sandwich ELISA of ab9667. Allows the detection of 0.2 - 0.4 ng/well of recombinant hGM-CSF (using 100µl/well antibody solution). Can be paired for Sandwich ELISA with [ab84259](#). Can be used as the capture antibody when paired with recommended pair.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors