

Anti-HA tag antibody [HA-7] ab49969

★★★★★ [5 Abreviews](#) [11 References](#) [画像数 3](#)

製品の概要

製品名	Anti-HA tag antibody [HA-7]
製品の詳細	Mouse monoclonal [HA-7] to HA tag
由来種	Mouse
アプリケーション	適用あり: ELISA, ICC/IF, IP, WB
種交差性	交差種: Species independent
免疫原	Synthetic peptide: YPYDVPDYA conjugated to KLH, corresponding to amino acids 98-106 of Influenza A Virus Hemagglutinin HA tag 📄 Run BLAST with 📄 Run BLAST with
ポジティブ・コントロール	Whole extracts of transfected 293T (human embryonal kidney) cells, expressing HA-tagged fusion protein.
特記事項	<p>This product was changed from ascites to tissue culture supernatant on 25th October 2016. The following lot(s) is/are from ascites and is still in stock as of 25th October 2016 - GR275412, GR177296. Lot numbers other than GR275412 and GR177296 will be from tissue culture supernatant. Please note that the dilutions may need to be adjusted accordingly.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.40 Preservative: 0.097% Sodium azide

	Constituents: 1% BSA, PBS
精製度	Proprietary Purification
特記事項 (精製)	Purified from Tissue culture supernatant".
ポリ/モノ	モノクローナル
クローン名	HA-7
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab49969の使用に適用されます**


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
ICC/IF		Use a concentration of 1 - 2 µg/ml.
IP	★★★★★ (1)	Use a concentration of 0.5 - 2 µg/ml.
WB	★★★★★ (2)	Use a concentration of 0.25 - 1 µg/ml. Predicted molecular weight: 63 kDa.

ターゲット情報


関連性	Human influenza hemagglutinin (HA) is a surface glycoprotein required for the infectivity of the human virus. The HA tag is derived from the HA molecule corresponding to amino acids 98-106 has been extensively used as a general epitope tag in expression vectors. Many recombinant proteins have been engineered to express the HA tag, which does not appear to interfere with the bioactivity or the biodistribution of the recombinant protein. This tag facilitates the detection, isolation, and purification of the proteins.
-----	--

画像



Immunocytochemistry/ Immunofluorescence - Anti-HA tag antibody [HA-7] (ab49969)

HEK-293 cells overexpressing N-terminal HA tagged fusion protein were fixed and permeabilized with 4% paraformaldehyde followed by 0.5% Triton™ X-100. Fixed cells were stained with 1 µg/mL Anti-HA tag antibody [HA-7] (ab49969). The antibody was developed using Goat Anti-Mouse IgG, Cy3™ conjugate. Cells were counterstained with DAPI (blue) to stain nuclei.




Immunoprecipitation - Anti-HA tag antibody [HA-7] (ab49969)

This image is courtesy of an anonymous Abreview

ab49969 Immunoprecipitating HA tag in HA tagged protein transfected 293T whole cell lysate. 25µg of cell lysate was incubated with primary antibody (1/2500 in lysis buffer) and matrix (Protein A/G) for 6 hours at 4°C. For western blotting **anti-mouse HRP (ab6789)** (1/10000) was used to confirm successful immunoprecipitation.

75kD is specific band for HA tagged protein. 120kD, 55kD, 35kD is non specific band for protein A/G bead.



Immunocytochemistry/ Immunofluorescence - Anti-HA tag antibody [HA-7] (ab49969)

HEK-293 cells overexpressing C-terminal HA tagged fusion protein were fixed and permeabilized with 4% paraformaldehyde followed by 0.5% Triton™ X-100. Fixed cells were stained with 1 µg/mL Anti-HA tag antibody [HA-7] (ab49969). The antibody was developed using Goat Anti-Mouse IgG, Cy3™ conjugate. Cells were counterstained with DAPI (blue) to stain nuclei.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors