

Anti-Hairy Cell Leukemia antibody [DBA.44], prediluted ab857

画像数 1

製品の概要

製品名	Anti-Hairy Cell Leukemia antibody [DBA.44], prediluted
製品の詳細	Mouse monoclonal [DBA.44] to Hairy Cell Leukemia, prediluted
由来種	Mouse
特異性	This antibody reacts with over 97% of hairy cell leukaemias as well as about 35% of high grade B cell lymphomas. It reacts with B lymphocytes, cells of mantle zone and immunoblasts outside the lymphoid tissues.
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	Tissue, cells or virus corresponding to Hairy Cell Leukemia. BALB/C mice were injected with the Deau cell line established from a large cell lymphoma.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Do Not Freeze.
バッファー	pH: 7.3 Preservative: 0.05% Sodium azide
ポリ/モノ	Inert stabilizer モノクローナル
クローン名	DBA.44
ミエローマ	unknown
アイソタイプ	IgM

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab857の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
IHC-P		Use at an assay dependent concentration.

ターゲット情報

細胞内局在

Nuclear

画像



Immunohistochemistry (Formalin fixed paraffin-embedded sections) analysis of human tonsil tissue labelling Hairy Cell Leukemia with ab857.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Hairy Cell Leukemia antibody [DBA.44], prediluted (ab857)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors