

Anti-Histone H2A.Z antibody [EPR6172] - ChIP Grade - BSA and Azide free ab248007

リコンビナント **RabMAb**

画像数 3

製品の概要

製品名	Anti-Histone H2A.Z antibody [EPR6172] - ChIP Grade - BSA and Azide free
製品の詳細	Rabbit monoclonal [EPR6172] to Histone H2A.Z - ChIP Grade – BSA and Azide free
由来種	Rabbit
特異性	This antibody is known to detect both isoforms of H2A.Z (H2A.Z.1 and H2A.Z.2).
アプリケーション	適用あり: WB, IHC-P, ChIP-sequencing 適用なし: Flow Cyt, ICC/IF or IP
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: C6, Raw 264.7, PC-12, NIH/3T3, HeLa, HepG2, and Human thymus lysates; IHC-P: Human breast tissue; ChIP-seq: Chromatin from HeLa cells.
特記事項	<p>ab248007 is the carrier-free version of ab124793.</p> <p>Our carrier-free antibodies are typically supplied in a PBS-only formulation, purified and free of BSA, sodium azide and glycerol. The carrier-free buffer and high concentration allow for increased conjugation efficiency.</p> <p>This conjugation-ready format is designed for use with fluorochromes, metal isotopes, oligonucleotides, and enzymes, which makes them ideal for antibody labelling, functional and cell-based assays, flow-based assays (e.g. mass cytometry) and Multiplex Imaging applications.</p> <p>Use our conjugation kits for antibody conjugates that are ready-to-use in as little as 20 minutes with <1 minute hands-on-time and 100% antibody recovery: available for fluorescent dyes, HRP, biotin and gold.</p> <p>This product is compatible with the Maxpar[®] Antibody Labeling Kit from Fluidigm, without the need for antibody preparation. Maxpar[®] is a trademark of Fluidigm Canada Inc.</p> <p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p>

Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to [RabMAb[®] patents](#).

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Do Not Freeze.
バッファー	pH: 7.2 Constituent: PBS
キャリア・フリー	はい
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR6172
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab248007の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent concentration. Predicted molecular weight: 14 kDa.
IHC-P		Use at an assay dependent concentration. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. Antigen retrieval is recommended.
ChIP-sequencing		Use at an assay dependent concentration.

追加情報 Is unsuitable for Flow Cyt, ICC/IF or IP.


ターゲット情報

機能	Variant histone H2A which replaces conventional H2A in a subset of nucleosomes. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability. DNA accessibility is regulated via a complex set of post-translational modifications of histones, also called histone code, and nucleosome remodeling. May be involved in the formation of constitutive heterochromatin. May be required for chromosome segregation during cell division.
配列類似性	Belongs to the histone H2A family.
翻訳後修飾	Monoubiquitination of Lys-122 gives a specific tag for epigenetic transcriptional repression.

Acetylated on Lys-5, Lys-8 and Lys-12 during interphase. Acetylation disappears at mitosis.
Monomethylated on Lys-5 and Lys-8 by SETD6. SETD6 predominantly methylates Lys-8, lys-5 being a possible secondary site.
Not phosphorylated.
Nucleus. Chromosome.

細胞内局在

画像



All lanes : Anti-Histone H2A.Z antibody [EPR6172] - ChIP Grade ([ab124793](#)) at 1/1000 dilution

Lane 1 : C6 lysates

Lane 2 : Raw 264.7 lysates

Lane 3 : PC-12 lysates

Lane 4 : NIH/3T3 lysates

Lane 5 : HeLa lysates

Lane 6 : HepG2 lysates

Lane 7 : Human thymus lysates


Lysates/proteins at 10 µg per lane.

Secondary

All lanes : HRP labelled goat anti-rabbit at 1/2000 dilution

Predicted band size: 14 kDa

This data was developed using [ab124793](#), the same antibody clone in a different buffer formulation.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Histone H2A.Z antibody [EPR6172] - ChIP Grade - BSA and Azide free (ab248007)

This data was developed using **ab124793**, the same antibody clone in a different buffer formulation. **ab124793**, at a 1/250 dilution, staining Histone H2A.Z in paraffin embedded Human breast tissue by Immunohistochemistry. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-Histone H2A.Z antibody [EPR6172] - ChIP Grade - BSA and Azide free (ab248007)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors