

Product datasheet

Anti-HPV16 L1 antibody [HPV16/2058R] ab234305

リコンビナント

画像数 2

製品の概要

製品名	Anti-HPV16 L1 antibody [HPV16/2058R]
製品の詳細	Rabbit monoclonal [HPV16/2058R] to HPV16 L1
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human papillomavirus
免疫原	Recombinant full length protein corresponding to Human papillomavirus HPV16 L1. Database link: P03101
ポジティブ・コントロール	IHC-P: Human cervix tissue.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 7.2 Preservative: 0.05% Sodium azide Constituents: PBS, 0.05% BSA
精製度	Protein A/G purified
特記事項(精製)	Purified from bioreactor concentrate by Protein A/G.
ポリ/モノ	モノクローナル
クローン名	HPV16/2058R
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab234305の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-P		Use a concentration of 0.5 - 1 µg/ml. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol. Incubate with primary antibody for 30 minutes at room temperature.

ターゲット情報

関連性	L1 is a major capsid protein of type 16 human papilloma virus. Infection with specific types of HPV has been associated with an increased risk of developing cervical neoplasia. HPV types 6 and 11 have been associated with relatively benign diseases such as genital warts but types 16 and 18 are strongly associated with cervical, vaginal, and vulvar malignancies.
細胞内局在	Virion. Host Nucleus.

画像



Formalin-fixed, paraffin-embedded human cervix tissue stained for HPV16 L1 using ab234305 at 1 µg/ml in immunohistochemical analysis.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-HPV16 L1 antibody [HPV16/2058R] (ab234305)

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-HPV16 L1 antibody [HPV16/2058R] (ab234305)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors