



HRP Anti-HA tag antibody ab183884

1 References [画像数 1](#)

製品の概要

| | |
|--------------|---|
| 製品名 | HRP Anti-HA tag antibody |
| 製品の詳細 | HRP Rabbit polyclonal to HA tag |
| 由来種 | Rabbit |
| 標識 | HRP |
| アプリケーション | 適用あり: WB |
| 免疫原 | Synthetic peptide corresponding to HA tag aa 100-200 conjugated to keyhole limpet haemocyanin. Database link: P03437 |
| |  Run BLAST with  Run BLAST with |
| ポジティブ・コントロール | Transfected HeLa cells overexpressing HA-tagged protein. |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|--------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot. Store at +4°C. Avoid freeze / thaw cycle. |
| バッファー | pH: 7.40 Preservative: 0.025% Proclin 300 Constituents: 1% BSA, 0.88% Sodium chloride, 0.6% Tris |
| 精製度 | Protein A purified |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab183884の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---------|
| WB | | 1/3000. |

ターゲット情報

関連性

Human influenza hemagglutinin (HA) is a surface glycoprotein required for the infectivity of the human virus. The HA tag is derived from the HA molecule corresponding to amino acids 98-106 has been extensively used as a general epitope tag in expression vectors. Many recombinant proteins have been engineered to express the HA tag, which does not appear to interfere with the bioactivity or the biodistribution of the recombinant protein. This tag facilitates the detection, isolation, and purification of the proteins.

画像



All lanes : HRP Anti-HA tag antibody (ab183884) at 1/3000 dilution

Lane 1 : Non-transfected HeLa cell lysate

Lane 2 : Transfected HeLa cells over expressing HA-tagged protein lysate

7.5% SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors