

HRP Anti-HLA Class I antibody [W6/32] ab33256

2 References

製品の概要

製品名	HRP Anti-HLA Class I antibody [W6/32]
製品の詳細	HRP Mouse monoclonal [W6/32] to HLA Class I
由来種	Mouse
標識	HRP
特異性	This antibody recognises an antigenic determinant shared among products of the HLA A, B and C loci.
アプリケーション	適用あり: IHC-Fr 適用なし: IHC-P
種交差性	交差種: Cow, Cat, Human, Monkey, Baboon, Cynomolgus monkey, Macaque monkey
免疫原	Tissue, cells or virus corresponding to HLA Class I. Purified human tonsil lymphocyte membranes.
エピトープ	Clone W6/32 recognises a conformational epitope, reacting with HLA class I alpha 3 and alpha 2 domains.
ポジティブ・コントロール	Tonsil.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	pH: 7.40 Preservative: 0.09% Sodium azide Constituent: PBS
精製度	Protein G purified
特記事項 (精製)	This antibody was purified from tissue culture supernatant.

ポリ/モノ	モノクローナル
クローン名	W6/32
ミエローマ	NS1/1-Ag4-1
アイソタイプ	IgG2a

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab33256の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IHC-Fr		Use at an assay dependent concentration.

追加情報 Is unsuitable for IHC-P.

ターゲット情報

関連性 HLA CClass I is involved in the presentation of foreign antigens to the immune system.
細胞内局在 Plasma membrane

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors