

HRP Anti-IL-7 antibody ab106025

画像数 1

製品の概要

製品名	HRP Anti-IL-7 antibody
製品の詳細	HRP Rabbit polyclonal to IL-7
由来種	Rabbit
標識	HRP
特異性	ab106025 will recognize both native and recombinant Human IL7 in cell supernatants and certain body fluids.
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant full length protein corresponding to Human IL-7.
ポジティブ・コントロール	Recombinant Human IL7
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.01% Gentamicin sulphate Constituents: 0.42% Potassium phosphate, 0.88% Sodium chloride, 1% BSA
精製度	Do NOT add Sodium Azide! Ion Exchange Chromatography
特記事項 (精製)	ab106025 is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer. This purified antibody has been heated to 56°C for 30 minutes.

ポリ/モノ
アイソタイプ

ポリクローナル
IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab106025の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/5000. Detects a band of approximately 17 kDa (predicted molecular weight: 20 kDa).


ターゲット情報

機能 Hematopoietic growth factor capable of stimulating the proliferation of lymphoid progenitors. It is important for proliferation during certain stages of B-cell maturation.

配列類似性 Belongs to the IL-7/IL-9 family.

細胞内局在 Secreted.

画像



HRP Anti-IL-7 antibody (ab106025) at 1/1000 dilution +
Recombinant Human IL7

Secondary

peroxidase conjugated anti-Rabbit IgG at 1/40000 dilution

Predicted band size: 20 kDa

Observed band size: 17 kDa

Additional bands at: 150 kDa (possible dimer)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors