

Anti-Hsp101 antibody - C-terminal ab80121

[2 References](#) [画像数 1](#)

製品の概要

製品名	Anti-Hsp101 antibody - C-terminal
製品の詳細	Rabbit polyclonal to Hsp101 - C-terminal
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Arabidopsis thaliana
免疫原	Recombinant fragment from the C-terminal of <i>Arabidopsis thaliana</i> Hsp101.
ポジティブ・コントロール	Arabidopsis thaliana WT lysate stressed at 38°C for 1.5 hours.
特記事項	<p>ab80121 also reacts with <i>Agave tequilana</i> (monocot).</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Lyophilized:For reconstitution add 200ul of sterile water.
保存方法	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
バッファー	Constituent: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab80121の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		1/1000. Predicted molecular weight: 101 kDa. 1/1000 with ECL.

ターゲット情報

関連性 Hsps/chaperones are known to be expressed in plants not only when they experience high temperatures stress, but also in response to a wide range of other environmental stresses such as water stress, salt stress, cold stress and oxidative stress. Hsp100 is one such major class that acts as a chaperone. Hsp101 is a member of this Hsp100 family. These proteins help dissociate protein aggregates formed during heat stress to allow them to be refolded by other chaperones. In spite of expression during heat stress, members of Hsp100 protein family are also expressed during seed development.

細胞内局在 Cytoplasmic

画像



Anti-Hsp101 antibody - C-terminal (ab80121) + Lysates prepared from Arabidopsis species.

Predicted band size: 101 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors