

Anti-IL-10 antibody [EP1115Y] ab52909

リコンビナント RabMAb[®]

3 References [画像数 2](#)

製品の概要

製品名	Anti-IL-10 antibody [EP1115Y]
製品の詳細	Rabbit monoclonal [EP1115Y] to IL-10
由来種	Rabbit
アプリケーション	適用あり: WB 適用なし: ICC/IF or IHC-P
種交差性	交差種: Human
免疫原	Synthetic peptide within Human IL-10 aa 100-200 (N terminal). The exact sequence is proprietary.
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none">- High batch-to-batch consistency and reproducibility- Improved sensitivity and specificity- Long-term security of supply- Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 0.31% Sodium citrate, 0.175% Sodium chloride, 0.0172% EDTA disodium salt, 59% PBS, 40% Glycerol (glycerin, glycerine), 0.05% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EP1115Y
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab52909の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		1/5000. Detects a band of approximately 21 kDa (predicted molecular weight: 21 kDa). This has only been tested against the recombinant protein.

追加情報 Is unsuitable for ICC/IF or IHC-P.

ターゲット情報

機能	Inhibits the synthesis of a number of cytokines, including IFN-gamma, IL-2, IL-3, TNF and GM-CSF produced by activated macrophages and by helper T-cells.
組織特異性	Produced by a variety of cell lines, including T-cells, macrophages, mast cells and other cell types.
配列類似性	Belongs to the IL-10 family.
細胞内局在	Secreted.

画像



Anti-IL-10 antibody [EP1115Y] (ab52909) at 1/5000 dilution +
Recombinant protein (5ng)

Secondary

Goat anti-Rabbit HRP labeled at 1/2000 dilution

Predicted band size: 21 kDa

Observed band size: 21 kDa

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-IL-10 antibody [EP1115Y] (ab52909)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors