

Anti-Insulin + Proinsulin antibody [D6C4] ab8304

★★★★★ **6 Abreviews** **30 References**

製品の概要

製品名	Anti-Insulin + Proinsulin antibody [D6C4]
製品の詳細	Mouse monoclonal [D6C4] to Insulin + Proinsulin
由来種	Mouse
特異性	Specific for insulin and proinsulin.
アプリケーション	適用あり: ELISA
種交差性	交差種: Mouse, Rat, Cow, Human, Pig
免疫原	Human insulin.
特記事項	<p>Kd for this antibody is 8.1×10^{-8}M.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	<p>pH: 7.40</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituent: 99.91% PBS</p>
精製度	Protein A purified
一次抗体 備考	Detection of both insulin and proinsulin.
ポリ/モノ	モノクローナル
クローン名	D6C4
ミエローマ	unknown
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab8304の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA	★★★★★ (1)	Use a concentration of 10 µg/ml.

ターゲット情報

関連性 Insulin decreases blood glucose concentration. It increases cell permeability to monosaccharides, amino acids and fatty acids. It accelerates glycolysis, the pentose phosphate cycle, and glycogen synthesis in liver. Defects in insulin are the cause of familial hyperproinsulinemia.

細胞内局在 Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors