

Anti-Integrin alpha 5 antibody [EPR7854] ab150361

KO 評価済 リコンビナント RabMAb

★★★★★ 12 Abreviews 68 References 画像数 18

製品の概要

製品名	Anti-Integrin alpha 5 antibody [EPR7854]
製品の詳細	Rabbit monoclonal [EPR7854] to Integrin alpha 5
由来種	Rabbit
アプリケーション	適用あり: Flow Cyt (Intra), WB, IHC-P, ICC/IF, IP
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide within Human Integrin alpha 5 aa 1000-1100 (C terminal). The exact sequence is proprietary.
ポジティブ・コントロール	WB: Human bladder, HT-1080, HeLa, U937, NIH/3T3, MEF, C2C12, Rat-1, and Rat liver lysates. IHC-P: Human kidney, human cerebral cortex, mouse cerebral cortex, and rat cerebral cortex tissues. ICC/IF: U937, MCF7 and wild-type HAP1 cells. Flow Cyt (intra): HeLa cells. IP: HeLa cells.
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb® technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb® patents.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.01% Sodium azide Constituents: 59% PBS, 40% Glycerol, 0.05% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR7854

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab150361の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Flow Cyt (Intra)		1/20. ab172730 - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody. For unpurified use at 1/100 - 1/100.
WB	★★★★★ (6)	1/1000 - 1/10000. Detects a band of approximately 150 kDa (predicted molecular weight: 115 kDa).
IHC-P	★★★★★ (2)	1/400. Perform heat mediated antigen retrieval before commencing with IHC staining protocol. For unpurified use at 1/100 - 1/250. See IHC antigen retrieval protocols
ICC/IF	★★★★★ (4)	1/50. For unpurified use at 1/250 - 1/500.
IP		1/10 - 1/100.

ターゲット情報

機能	Integrin alpha-5/beta-1 is a receptor for fibronectin and fibrinogen. It recognizes the sequence R-G-D in its ligands. In case of HIV-1 infection, the interaction with extracellular viral Tat protein seems to enhance angiogenesis in Kaposi's sarcoma lesions.
配列類似性	Belongs to the integrin alpha chain family. Contains 7 FG-GAP repeats.
細胞内局在	Membrane.

画像



All lanes : Anti-Integrin alpha 5 antibody [EPR7854] (ab150361) at 1/2000 dilution (Unpurified)

Lane 1 : Human bladder tissue lysate

Lane 2 : HT-1080 (Human urinary bladder carcinoma cell line) whole cell lysates

Lane 3 : HeLa (Human epithelial cell line from cervix adenocarcinoma) whole cell lysates

Lane 4 : U937 (Human histiocytic lymphoma cell line) whole cell lysates

Lysates/proteins at 10 µg per lane.

Secondary


All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size: 115 kDa

Observed band size : 19, 150 kDa


Exposure time : Lane 1-2: 1 minutes, Lane 3-4: 3 minutes

Blocking/Diluting buffer 5% NFDM/TBST



Immunocytochemistry/ Immunofluorescence - Anti-Integrin alpha 5 antibody [EPR7854] (ab150361)

Immunocytochemistry/ Immunofluorescence analysis of U937 (Human histiocytic lymphoma monocyte) cells labeling Integrin alpha 5 with purified ab150361 at 1:50 dilution (4.3 µg/ml). Cells were fixed in 4% Paraformaldehyde and permeabilized with 0.1% tritonX-100. Cells were counterstained with Ab195889 Anti-alpha Tubulin antibody [DM1A] - Microtubule Marker (Alexa Fluor® 594) 1:200 (2.5 µg/ml). Goat anti rabbit IgG (Alexa Fluor® 488, **ab150077**) was used as the secondary antibody at 1:1000 (2 µg/ml) dilution. DAPI (blue) was used as nuclear counterstain. PBS instead of the primary antibody was used as the secondary antibody only control.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Integrin alpha 5 antibody [EPR7854] (ab150361)

Immunohistochemical analysis of paraffin embedded Human kidney tissue labelling Integrin alpha 5 with unpurified ab150361 antibody at a dilution of 1/100.

Perform heat mediated antigen retrieval before commencing with IHC staining protocol.



Immunoprecipitation - Anti-Integrin alpha 5 antibody
[EPR7854] (ab150361)

ab150361 (Purified) at 1:20 dilution (1 µg) immunoprecipitating Integrin alpha 5 in HeLa whole cell lysate.


Lane 1 (input): HeLa (Human epithelial cell line from cervix adenocarcinoma) whole cell lysate 10 µg

Lane 2 (+): ab150361 & HeLa whole cell lysate

Lane 3 (-): Rabbit monoclonal IgG (**ab172730**) instead of ab150361 in HeLa whole cell lysate


For western blotting, VeriBlot for IP Detection Reagent (HRP) (**ab131366**) was used at 1:5000 dilution.

Blocking and diluting buffer: 5% NFDM/TBST.



Flow Cytometry (Intracellular) - Anti-Integrin alpha 5
antibody [EPR7854] (ab150361)

Intracellular Flow Cytometry analysis of HeLa (Human epithelial cell line from cervix adenocarcinoma) cells labeling Integrin alpha 5 with purified ab150361 at 1/20 dilution (10 µg/ml) (Red). Cells were fixed with 4% Paraformaldehyde and permeabilised with 90% Methanol. A Goat anti rabbit IgG (Alexa Fluor® 488, **ab150077**) secondary antibody was used at 1/2000. Isotype control - Rabbit monoclonal IgG (Black). Unlabeled control - Cell without incubation with primary antibody and secondary antibody (Blue).



All lanes : Anti-Integrin alpha 5 antibody [EPR7854] (ab150361) at 1/1000 dilution (Purified)

Lane 1 : NIH/3T3 (Mouse embryo fibroblast cell line) whole cell lysates

Lane 2 : MEF (Mouse embryonic fibroblast cell line) whole cell lysates

Lane 3 : C2C12 (Mouse myoblast cell line) whole cell lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution


Predicted band size: 115 kDa

Observed band size: 150,19 kDa

The bands detected by ab150361 are consistent with what have been described in PMID: 22740495. Full length: 150kDa; Light chain: 19kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Integrin alpha 5 antibody [EPR7854] (ab150361)



Western blot - Anti-Integrin alpha 5 antibody [EPR7854] (ab150361)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of mouse cerebral cortex tissue sections labeling Integrin alpha 5 with purified ab150361 at 1:400 dilution (0.54 µg/ml). Heat mediated antigen retrieval was performed using **ab93684** (Tris/EDTA buffer, pH 9.0). ImmunoHistoProbe one step HRP Polymer (ready to use) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain.

All lanes : Anti-Integrin alpha 5 antibody [EPR7854] (ab150361) at 1/1000 dilution (Purified)

Lane 1 : Rat1 (Rat fibroblast cell line) whole cell lysates

Lane 2 : Rat liver lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution

Predicted band size: 115 kDa


Observed band size: 150,19 kDa

The bands detected by ab150361 are consistent with what have been described in PMID: 22740495. Full length: 150kDa; Light chain: 19kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Integrin alpha 5 antibody [EPR7854] (ab150361)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of rat cerebral cortex tissue sections labeling Integrin alpha 5 with purified ab150361 at 1:400 dilution (0.54 µg/ml). Heat mediated antigen retrieval was performed using **ab93684** (Tris/EDTA buffer, pH 9.0). ImmunoHistoProbe one step HRP Polymer (ready to use) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain.



Western blot - Anti-Integrin alpha 5 antibody [EPR7854] (ab150361)

Anti-Integrin alpha 5 antibody [EPR7854] (ab150361) at 1/10000 dilution (Purified) + HT-1080 (Human fibrosarcoma cell line) whole cell lysates at 15 µg


Secondary

Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution

Predicted band size: 115 kDa

Observed band size: 150,19 kDa

The bands detected by ab150361 are consistent with what have been described in PMID: 22740495. Full length: 150kDa; Light chain: 19kDa



Western blot - Anti-Integrin alpha 5 antibody
[EPR7854] (ab150361)

Anti-Integrin alpha 5 antibody [EPR7854] (ab150361) at 1/1000 dilution (Purified) + U937 (Human histiocytic lymphoma cell line) whole cell lysates at 15 µg


Secondary

Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution

Predicted band size: 115 kDa


Observed band size: 150,19 kDa

The bands detected by ab150361 are consistent with what have been described in PMID: 22740495. Full length: 150kDa; Light chain: 19kDa




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Integrin alpha 5 antibody
[EPR7854] (ab150361)

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) analysis of human cerebral cortex tissue sections labeling Integrin alpha 5 with purified ab150361 at 1:400 dilution (0.54 µg/ml). Heat mediated antigen retrieval was performed using [ab93684](#) (Tris/EDTA buffer, pH 9.0). ImmunoHistoProbe one step HRP Polymer (ready to use) was used as the secondary antibody. Negative control: PBS instead of the primary antibody. Hematoxylin was used as a counterstain.




Immunocytochemistry/ Immunofluorescence - Anti-
Integrin alpha 5 antibody [EPR7854] (ab150361)

ab150361 staining Integrin alpha 5 in MCF7 (Human breast adenocarcinoma cell line) cells. The cells were fixed with 100% methanol (5min), permeabilized with 0.1% Tween for 5 minutes and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Tween for 1h. The cells were then incubated with unpurified ab150361 at 1 µg/ml concentration and **ab7291** at 1 µg/ml concentration overnight at +4°C, followed by a further incubation at room temperature for 1h with a goat secondary antibody to Rabbit IgG (Alexa Fluor® 488) (**ab150081**) at 2 µg/ml (shown in green) and a goat secondary antibody to Mouse IgG (Alexa Fluor® 594) (**ab150120**) at 2 µg/ml (showed in pseudocolour red) . Nuclear DNA was labelled in blue with DAPI.




Immunocytochemistry/ Immunofluorescence - Anti-
Integrin alpha 5 antibody [EPR7854] (ab150361)

unpurified ab150361 staining Integrin alpha 5 in wild-type HAP1 cells (top panel) and ITGA5 knockout HAP1 cells (bottom panel). The cells were fixed with 100% methanol (5min), permeabilized with 0.1% Tween for 5 minutes and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Tween for 1h. The cells were then incubated with ab150361 at 1/250 dilution and **ab7291** at 1 µg/ml concentration overnight at +4°C, followed by a further incubation at room temperature for 1h with a goat secondary antibody to Rabbit IgG (Alexa Fluor® 488) (**ab150081**) at 2 µg/ml (shown in green) and a goat secondary antibody to Mouse IgG (Alexa Fluor® 594) (**ab150120**) at 2 µg/ml (showed in pseudocolour red) . Nuclear DNA was labelled in blue with DAPI.




Immunocytochemistry/ Immunofluorescence - Anti-Integrin alpha 5 antibody [EPR7854] (ab150361)

Immunofluorescent analysis of 4% paraformaldehyde-fixed, 0.1% Triton X-100 permeabilized U937 (Human histiocytic lymphoma cell line) cells labeling Integrin alpha 5 with unpurified ab150361 at 1/500 dilution, followed by Goat anti-rabbit IgG (Alexa Fluor 488) (**ab150077**) secondary antibody at 1/400 dilution (green). Confocal image showing mainly membrane staining on U937 cell line. The nuclear counterstain is DAPI (blue).



Flow Cytometry (Intracellular) - Anti-Integrin alpha 5 antibody [EPR7854] (ab150361)

Overlay histogram showing HeLa cells stained with unpurified ab150361 (red line). The cells were fixed with 4% paraformaldehyde (10 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab150361, 1/1000 dilution) for 30 min at 22°C. The secondary antibody used was goat anti-rabbit Alexa Fluor 488 IgG (H&L) (**ab150077**) at 1/2000 dilution for 30 min at 22°C. Isotype control antibody (black line) was rabbit IgG (monoclonal) (0.1 µg/1x10⁶ cells) used under the same conditions. Unlabelled sample (blue line) was also used as a control. Acquisition of >5,000 events were collected using a 20mW Argon ion laser (488nm) and 525/30 bandpass filter. This antibody gave a positive signal in HeLa cells fixed with 80% methanol (5 min)/permeabilized with 0.1% PBS-Tween for 20 min used under the same conditions.



Western blot - Anti-Integrin alpha 5 antibody [EPR7854] (ab150361)

All lanes : Anti-Integrin alpha 5 antibody [EPR7854] (ab150361) at 1/2000 dilution (Unpurified)

Lane 1 : Mouse heart tissue lysate

Lane 2 : Rat heart tissue lysate

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/1000 dilution

Predicted band size: 115 kDa

Exposure time: 3 minutes

Observed band size : 19, 150 kDa

Exposure time : 3 minutes

Blocking/Diluting buffer 5% NFDm/TBST

Why choose a recombinant antibody?



Anti-Integrin alpha 5 antibody [EPR7854] (ab150361)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors