

# Anti-Interferon alpha 1 antibody [ST29] ab20200

**2 References**   [画像数 1](#)

### 製品の概要

製品名	Anti-Interferon alpha 1 antibody [ST29]
製品の詳細	Mouse monoclonal [ST29] to Interferon alpha 1
由来種	Mouse
アプリケーション	<b>適用あり:</b> Sandwich ELISA
種交差性	<b>交差種:</b> Recombinant fragment
免疫原	Human interferon (HuIFN) produced by Namalwa cells (HuIFN-aN).
特記事項	Myeloma: NS2  The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.  If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
精製度	Protein A/G purified
ポリ/モノ	モノクローナル
クローン名	ST29
アイソタイプ	IgG1
軽鎖の種類	unknown

### アプリケーション


**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab20200の使用に適用されます  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
Sandwich ELISA		Use a concentration of 5 µg/ml.

## ターゲット情報

機能	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.
配列類似性	Belongs to the alpha/beta interferon family.
細胞内局在	Secreted.

## 画像



Standard Curve for Interferon alpha (Analyte: **ab48750**) dilution range 1pg/ml to 1ug/ml using Capture Antibody Mouse monoclonal [ST29] to Interferon alpha (ab20200) at 5ug/ml and Detector Antibody Rabbit polyclonal to Interferon alpha (Biotin) (**ab84257**) at 0.5ug/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors