

Anti-Isopeptide antibody [81D1C2] ab422

★★★★★ **1 Abreviews** **10 References**

製品の概要

製品名	Anti-Isopeptide antibody [81D1C2]
製品の詳細	Mouse monoclonal [81D1C2] to Isopeptide
由来種	Mouse
特異性	This antibody reacts with the isopeptide : epsilon-(gamma-glutamyl) lysine. There is no cross reaction with either lysine or glutamine. They are specific to the epsilon-(gamma-glutamyl) lysine bond formed by transglutaminase activity.
アプリケーション	適用あり: WB
種交差性	交差種: Species independent
免疫原	Synthetic isopeptide.
特記事項	<p>The ability of this antibody to fix complement is unknown. The isopeptide may be purchased from Bachem Catalog number G-1970.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Constituents: 0.75% Glycine, 1.21% Tris, 2% Sucrose
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	81D1C2
ミエローマ	unknown
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab422の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		1/50 - 1/500.

ターゲット情報

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors