

Anti-LC3B antibody ab51520

★★★★★ **18 Abreviews** **399 References** 画像数 9

製品の概要

製品名	Anti-LC3B antibody
製品の詳細	Rabbit polyclonal to LC3B
由来種	Rabbit
特異性	ab51520 is specific for the B isoform of LC3.
アプリケーション	適用あり: IHC-Fr, Flow Cyt (Intra), WB, ICC/IF
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to Human LC3B (N terminal). Synthetic peptide made to the N terminal region of human LC3, isoform B protein. Database link: Q9GZQ8
ポジティブ・コントロール	ICC/IF: HeLa, U373-MGc cells. IHC-P: U87-MG cells. Flow Cyt (Intra): NTERA-2 cells. IHC-P-Fr: U373-MG cells.
特記事項	<p>ab51520 has been shown to work on human malignant brain tumors, breast cancer, leukemia samples, and murine cells.</p> <p>The product formulation changed to include 1% BSA on 17th May 2017. The following lot does not contain BSA and is still in stock GR320876. All future lots will contain 1% BSA.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: PBS, 1% BSA
精製度	Immunogen affinity purified

一次抗体 備考	ab51520 has been shown to work on human malignant brain tumors, breast cancer, leukemia samples, and murine cells.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab51520の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
IHC-Fr	★★★★★ (2)	1/2000.
Flow Cyt (Intra)		1/200.
WB	★★★★★ (9)	1/3000. Predicted molecular weight: 15 kDa. Detects bands of ~16 kDa (LC3-II) and ~18 kDa (LC3-I).
ICC/IF	★★★★★ (3)	1/2000.

ターゲット情報

機能	Probably involved in formation of autophagosomal vacuoles (autophagosomes).
組織特異性	Most abundant in heart, brain, skeletal muscle and testis. Little expression observed in liver.
配列類似性	Belongs to the MAP1 LC3 family.
翻訳後修飾	The precursor molecule is cleaved by APG4B/ATG4B to form LC3-I. This is activated by APG7L/ATG7, transferred to ATG3 and conjugated to phospholipid to form LC3-II.
細胞内局在	Cytoplasm > cytoskeleton. Endomembrane system. Cytoplasmic vesicle > autophagosome membrane. LC3-II binds to the autophagic membranes.

画像



All lanes : Anti-LC3B antibody (ab51520) at 1/3000 dilution

Lane 1 : Wild-type HepG2 untreated control cell lysate

Lane 2 : Wild-type HepG2 Treated Chloroquine (50 uM, 16 h) cell lysate

Lane 3 : MAP1LC3B knockout HepG2 untreated control cell lysate

Lane 4 : MAP1LC3B knockout HepG2 Treated Chloroquine (50 uM, 16 h) cell lysate

Lane 5 : U-87 MG cell lysate

Lane 6 : PC-12 cell lysate


Lysates/proteins at 20 µg per lane.

Performed under reducing conditions.

Predicted band size: 15 kDa

Observed band size: 14,16 kDa


False colour image of Western blot: Anti-LC3B antibody staining at 1/3000 dilution, shown in green; Mouse anti-Alpha Tubulin [DM1A] ([ab7291](#)) loading control staining at 1/20000 dilution, shown in red. In Western blot, ab51520 was shown to bind specifically to LC3B. A band was observed at 16/14 kDa (yellow arrows) in treated wild-type HepG2 cell lysates with no signal observed at this size in MAP1LC3B knockout cell line [ab277828](#) (knockout cell lysate [ab283796](#)). To generate this image, wild-type and MAP1LC3B knockout HepG2 cell lysates were analysed. First, samples were run on an SDS-PAGE gel then transferred onto a nitrocellulose membrane. Membranes were blocked in 3 % milk in TBS-0.1 % Tween® 20 (TBS-T) before incubation with primary antibodies overnight at 4 °C. Blots were washed four times in TBS-T, incubated with secondary antibodies for 1 h at room temperature, washed again four times then imaged. Secondary antibodies used were Goat anti-Rabbit IgG H&L 800CW and Goat anti-Mouse IgG H&L 680RD at 1/20000 dilution.



Immunohistochemistry (Frozen sections) - Anti-LC3B antibody (ab51520)

Ab51520 at 1/2000 dilution staining treated U87-MG cultured and subcutaneous tumors

- 1) Control
- 2) Temozolomide



Western blot - Anti-LC3B antibody (ab51520)

All lanes : Anti-LC3B antibody (ab51520) at 1/3000 dilution


Lane 1 : treated U87-MG (human glioblastoma astrocytoma) lysates with Arsenic Trioxide at 0 μ M

Lane 2 : treated U87-MG (human glioblastoma astrocytoma) lysates with Arsenic Trioxide at 2 μ M

Lane 3 : treated U87-MG (human glioblastoma astrocytoma) lysates with Arsenic Trioxide at 4 μ M


Lane 4 : treated U87-MG (human glioblastoma astrocytoma) lysates with Arsenic Trioxide at 8 μ M

Predicted band size: 15 kDa




Immunocytochemistry/ Immunofluorescence - Anti-LC3B antibody (ab51520)

Immunocytochemistry/Immunofluorescence analysis of HeLa cells labelling LC3B with ab51520. An Alexa Fluor 488-conjugated Goat to rabbit IgG was used as secondary antibody (green). Actin filaments were labeled with Alexa Fluor 568 phalloidin (red). DAPI was used to stain the cell nuclei (blue).



Flow Cytometry (Intracellular) - Anti-LC3B antibody (ab51520)

Flow cytometry (Intracellular) analysis of NTERA-2 cells with ab51520 at 1/200 using Dylight-488 conjugated goat anti-rabbit IgG secondary.




Immunocytochemistry/ Immunofluorescence - Anti-LC3B antibody (ab51520)

Ab51520 at 1/2000 dilution staining treated U373-MG cells;

- 1) Control
- 2) C₂-ceramide.


The nuclei were stained with DAPI.



Immunohistochemistry (Frozen sections) - Anti-LC3B antibody (ab51520)


Ab51520 at 1/2000 dilution staining treated U373-MG (human glioblastoma) cells

- 1) Control
- 2) Temozolomide



Immunocytochemistry/ Immunofluorescence - Anti-LC3B antibody (ab51520)

Immunocytochemistry/Immunofluorescence analysis of HeLa cells using anti-LC3B antibody(red) at 0.1 ug/ml. Nuclei were counterstained with DAPI (blue).Nuclei and alpha-tubulin were counterstained with DAPI (blue) and Dylight 550 (red).



Immunocytochemistry/Immunofluorescence analysis of HeLa cells using anti-LC3B antibody(red). Nuclei were counterstained with DAPI (blue).

Immunocytochemistry/ Immunofluorescence - Anti-LC3B antibody (ab51520)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors