

Anti-Legionella pneumophila antibody ab20943

★★★★★ [1 Abreviews](#) [4 References](#)

製品の概要

製品名	Anti-Legionella pneumophila antibody
製品の詳細	Rabbit polyclonal to Legionella pneumophila
由来種	Rabbit
特異性	Reacts with L. pneumophila serogroups 1-12 in IFA. Recognizes all antigens of intact microorganism. Antiserum is not absorbed and may cross-react with related microorganisms
アプリケーション	適用あり: ELISA, ICC/IF, Conjugation
種交差性	交差種: Legionella pneumophila
免疫原	Tissue, cells or virus corresponding to Legionella pneumophila. A whole cell preparation of Legionella pneumophila; ATCC #33152.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab20943の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent concentration.
ICC/IF	★★★★★ (1)	Use at an assay dependent concentration.
Conjugation		Use at an assay dependent concentration.

ターゲット情報

関連性 Legionella pneumophila is a flagellated gram negative bacterium found primarily in warm water environments. Legionella pneumophila causes a type of pneumonia called Legionnaire disease and a milder condition called Pontiac fever. Infection is acquired through inhalation of the microorganism in water mists (e.g. from air conditioning, cooling towers, whirlpool spas and showers mists). A number of risk factors for acquiring Legionnaire disease have been identified, including age, smoking, chronic lung disease, cancer, and immunosuppression.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors