


### Anti-Ly-6A/E (Sca-1) antibody ab95439

★★★★☆ 1 Abreviews 2 References 画像数 1

**医薬用外劇物**

#### 製品の概要

製品名	Anti-Ly-6A/E (Sca-1) antibody
製品の詳細	Rabbit polyclonal to Ly-6A/E (Sca-1)
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Rat <b>交差が予測される動物種:</b> Mouse 
免疫原	Synthetic peptide, corresponding to amino acids 91-108 of Rat Sca1/ Ly6A/E
ポジティブ・コントロール	Rat heart, liver, testicular and brain tissue lysates.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservatives: 0.025% Sodium azide, 0.025% Thimerosal (merthiolate) Constituents: 2.5% BSA, 0.45% Sodium chloride, 0.1% Dibasic monohydrogen sodium phosphate
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

## The Abpromise guarantee

**Abpromise保証は、次のテスト済みアプリケーションにおけるab95439の使用に適用されず**


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 14 kDa.

## ターゲット情報

機能	T-cell activation.
組織特異性	Widely expressed.
配列類似性	Contains 1 UPAR/Ly6 domain.
翻訳後修飾	O-glycosylated. Not N-glycosylated. Not phosphorylated.
細胞内局在	Cell membrane.

## 画像



Western blot - Anti-Ly-6A/E (Sca-1) antibody (ab95439)

**All lanes :** Anti-Ly-6A/E (Sca-1) antibody (ab95439) at 1 µg/ml

**Lane 1 :** Rat heart tissue lysate

**Lane 2 :** Rat liver tissue lysate

**Lane 3 :** Rat testicular tissue lysate

**Lane 4 :** Rat brain tissue lysate

**Predicted band size:** 14 kDa

**Observed band size:** 16 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors