

Anti-mono methyl Arginine antibody [16B11] ab414

14 References

製品の概要

製品名	Anti-mono methyl Arginine antibody [16B11]
製品の詳細	Mouse monoclonal [16B11] to mono methyl Arginine
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Species independent
特記事項	<p>This antibody will be of central importance in analysing the methylation status of chromatin and transcription factors. Since it recognises monomethyl arginine only it can be used in parallel with ab413 (detects dimethyl arginine only) and ab412 (detects both monomethyl and dimethyl arginine) to monitor the exact modification status.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Constituents: 1.21% Tris, 0.75% Glycine, 2% Sucrose
精製度	Protein A purified
一次抗体 備考	<p>This antibody will be of central importance in analysing the methylation status of chromatin and transcription factors. Since it recognises monomethyl arginine only it can be used in parallel with ab413 (detects dimethyl arginine only) and ab412 (detects both monomethyl and dimethyl arginine) to monitor the exact modification status.</p>
ポリ/モノ	モノクローナル
クローン名	16B11

ミエローマ	unknown
アイソタイプ	IgG2b
軽鎖の種類	kappa

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab414の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/50 - 1/500.

ターゲット情報

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors